

Application note

NHS ester chemistry labeling protocol of Nanobodies

Content

1. Introduction.....	1
2. General considerations and recommendations.....	1
3. Preparation of dye	1
4. Preparation of VHH	2
5. Conjugation reaction.....	2
6. Removal of unbound dye.....	2
7. Storage	2

1. Introduction

This protocol provides recommendations for the random labeling of ChromoTek Nanobodies containing surface-exposed lysines with NHS-reactive fluorescent dyes by NHS ester chemistry. Note that Nanobodies containing ectopic cysteines are usually not suitable for NHS labeling.

2. General considerations and recommendations

- ▶ Each fluorescent dye is different and can influence the Nanobody to a different extent. The conditions for labeling must be established individually for each dye.
- ▶ Remember that Nanobodies are only 1/10 the size of an antibody when antibody labeling kits are used.
- ▶ Many fluorescent dyes have a hydrophobic structure. The conjugation of hydrophobic dyes to Nanobodies can affect the solubility of the Nanobody.

3. Preparation of dye

- ▶ Follow the dye manufacturer's protocol.
- ▶ Freshly prepare the dye stock solution immediately before starting the labeling reaction. Functional groups lose their reactivity during storage.

- ▶ Adjust the molar excess of the dye according to the dye manufacturer's recommendations. Aim for a degree of labeling of 1.
- ▶ Dyes are dissolved in organic solvents. Note that organic solvents can affect the stability and can facilitate precipitation of the Nanobody.

4. Preparation of VHH

- ▶ Centrifuge material before use (20,000x g, 15 min, +4°C).
- ▶ Nanobodies are stored in TAPS buffer with sodium azide (25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, 0.09 % sodium azide). Perform a buffer exchange step to TAPS buffer (25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA) to remove sodium azide. Note that the labeling buffer can influence the labeling efficiency.

5. Conjugation reaction

- ▶ Mix the diluted dye with the Nanobody.
- ▶ Incubate by end-over-end rotation at room temperature for 1-2 h.

6. Removal of unbound dye

- ▶ Centrifuge the solution after the labeling reaction is completed (20,000x g, 15 min, +4°C) and continue working with the supernatant.
- ▶ Separate unbound dye from the labelled Nanobody by one of the following options or by a combination thereof:
 - Size exclusion column (length: >30 cm)
 - Dialysis (molecular weight cut off: 3.5 kDa)
 - Spin column (molecular weight cut off: 7 kDa)
 - Desalting column

7. Storage

- ▶ Aliquot the labelled Nanobody and store at +4°C or -20°C. Avoid freeze-thaw cycles. Protect from light.
- ▶ Add 0.1% sodium azide for long-term storage to prevent bacterial contamination.

Disclaimer

ChromoTek is a registered trademark of ChromoTek GmbH, part of Proteintech Group. Nanobody is a registered trademark of Ablynx, a Sanofi company. Other suppliers' products may be trademarks or registered trademarks of the corresponding supplier each. Statements on other suppliers' products are given according to our best knowledge.