



EXPERTS IN SIZE AND CRYSTALS

Copenhagen Nanosystems ApS

TESTING AND VALIDATION OF NANOCUVETTE™ S, SPECTROWORKS™ TOGETHER WITH UV-VIS SPECTROPHOTOMETER

30-Jul-2021

WP



PARTICLE
ANALYTICAL

EXPERTS IN SIZE AND CRYSTALS

TESTING AND VALIDAITON OF NANOCUVETTE™ S, SPECTROWORKS™
TOGETHER WITH UV-VIS SPECTROPHOTOMETER

Date of Issue: 30-Jul-2021

Supersedes: -

Written by:

Wenbo Wang

Wenbo Wang, QC

Date:

30 JULI 2021

Approved by:

Søren Lund Kristensen

Søren Lund Kristensen, QA

Date:

30 JULI 2021

Sponsor:

Sponsor:

Name:

Name:

Date:

Date:

Sponsor:

Name:

Date:



TABLE OF CONTENTS

Abstract	4
1 Aim	4
2 Materials	4
3 Analytical equipment	5
4 Method background	5
4.1 UV-Vis measurements with NanoCuvette™ S and SpectroWorks™ modeling	5
4.2 DLS measurement	6
5 Analysis.....	6
5.1 Comparison of the diluted polystyrene latex beads concentration from stock provided by distributor with the concentration obtained from NanoCuvette™ S and SpectroWorks™	6
5.2 Comparison of the mean particle size of polystyrene latex beads from DLS and NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer...7	7
6 Conclusion	8

Appendix 1. Results from NanoCuvette™ S and SpectroWork™ together with UV-Vis spectrophotometer
(24 Pages)

Appendix 2. Results from dynamic light scattering (60 Pages)

wp 9



Abstract

The present report evaluates the particle size distribution and concentration of different polystyrene latex beads using NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer. The report includes:

- Evaluate mean particle size and concentration of polystyrene latex beads analyzed in SpectroWorks™ by using UV-Vis measurements with NanoCuvette™ S.
- Compare the result of concentration of polystyrene from NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer and diluted sample from stock provided by suppliers. Samples were prepared by diluting the stock provided by supplier (Sigma Aldrich).
- Compare the result of mean particle size of polystyrene latex beads from NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer and dynamic light scattering (DLS).

The performance of Nanocuvette™ S and Spectroworks™ is assessed by comparing the results of particle size and concentration of polystyrene latex beads done with DLS. Optical filter inserted cuvette, Nanocuvette™ S and software together with spectrophotometer provided reliable results where the concentrations of all samples (except for polystyrene latex beads 1100 nm) have deviation less than 30% compared to the concentration of sample dilutions. Furthermore, the particle size results obtained from analysis of the UV-Vis measurements with NanoCuvette™ S by SpectroWorks™ have less than 22% deviation (except for run 2 of 100 nm sample) from the standard for all the size range.

1 Aim

The present report documents the results of mean particle size and concentration of polystyrene latex beads analyzed by SpectroWorks™ using UV-Vis measurements with NanoCuvette™ S. To evaluate the reliability of the results obtained from the cphnano products, the results will be compared with DLS results and product specification obtained from sample dilutions and product specification sheets provided by suppliers (Sigma).

2 Materials

The material used for measurement with NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer are listed in Table 1.

Table 1: Sample received for UV-Vis measurement with NanoCuvette™ S and SpectroWorks™ modeling

Sample name	Concentration (Vol.) (diluted samples)
Polystyrene Latex beads 100 nm	0.01%
Polystyrene Latex beads 460 nm	0.003%
Polystyrene Latex beads 600 nm	0.003%
Polystyrene Latex beads 800 nm	0.003%
Polystyrene Latex beads 1100nm	0.001%
Polystyrene Latex beads 3000 nm	0.001%



The material used for DLS measurement are listed in Table 2.

Table 2: Samples received for testing and validation of NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer

Sample name	Batch number	Internal number
Polystyrene Latex beads 100 nm	2021/07/23	32027
Polystyrene Latex beads 460 nm	2021/07/23	32028
Polystyrene Latex beads 600 nm	2021/07/23	32029
Polystyrene Latex beads 800 nm	2021/07/23	32030
Polystyrene Latex beads 1100nm	2021/07/23	32031
Polystyrene Latex beads 3000 nm	2021/07/23	32032

3 Analytical equipment

Analytical equipment is listed in 3.

Table 3. Analytical equipment included in the project

UV-Vis spectrophotometer	VWR UV-6300 Double beam spectrophotometer
DLS unit	Malvern Zetasizer Nano ZS

4 Method background

4.1 UV-Vis measurements with NanoCuvette™ S and SpectroWorks™ modeling

UV-Vis measurements using NanoCuvette™ S was done from 200nm-1100nm. Firstly, reference measurement was done by pipetting 100 µL of reference solution, DI water in the NanoCuvette™ S. B, A and D side measurements (see figure 1) were done respectively, and the spectra were saved. Secondly, the reference sample was removed and 100 µL of diluted samples were added in the NanoCuvette™ S and B, A and D side spectra are recorded and saved for further analysis. Each measurement was replicated. Reference spectra A, B, D sides and sample spectra A, B and D sides were respectively dragged and dropped in the SpectroWorks™ to obtain the results.



4.2 DLS measurement

The measurement conditions for dynamic light scattering are described in Table 4.

Table 4. Measurement conditions of samples

Sample cell	Disposable standard 1 cm cuvette, volume 3 ml
Sample preparation	The sample is diluted Ca. 50 times: 20 uL of sample is added to 1 mL of water.
Dispersant refractive index	1.330
Viscosity	0.8872 cP
Material refractive index	1.590
Material Absorption	0.010
Temperature	25.0 °C
Equilibration time	180 seconds
Evaluation model	General purpose
Display range	0.600 to 6000 nm
Multimodal – analysis resolution	Normal
Lower threshold	0.05
Upper threshold	0.01
Measurement position	Automatic
Attenuator	Automatic
Run duration	Automatic
Number of runs per measurement	Automatic
Measurements per sampling	3
Delay between measurements	0 seconds
Number of samplings	2

5 Analysis

5.1 Comparison of the diluted polystyrene latex beads concentration from stock provided by distributor with the concentration obtained from NanoCuvette™ S and SpectroWorks™

The results obtained from analyzed data using SpectroWorks™ are listed in table 5





Table 5. Concentration of the diluted polystyrene latex beads obtained using NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer.

Sample name	Concentration (Vol.) Diluted from the stock	Concentration (Vol.) From NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer. Result (% deviation)	
		Run 1	Run 2
Polystyrene Latex beads 100 nm	0.01%	0.0099% (1.0%)	0.0081% (19.0%)
Polystyrene Latex beads 460 nm	0.003%	0.0038% (26.7%)	0.0038% (26.7%)
Polystyrene Latex beads 600 nm	0.003%	0.0026% (13.3%)	0.0026% (13.3%)
Polystyrene Latex beads 800 nm	0.003%	0.0030% (0.0%)	0.0031% (3.3%)
Polystyrene Latex beads 1100nm	0.001%	0.0016% (60.0%)	0.0016% (60.0%)
Polystyrene Latex beads 3000 nm	0.001%	0.0008% (20.0%)	0.0007% (30.0%)

The results of concentration obtained using NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer are considered reliable. The concentrations of all samples (except for polystyrene latex beads 1100 nm) from SpectroWroks™ have deviation less than 30%, when compared to the concentration of the dilution made from stock. However, due to the low concentration of polystyrene latex beads 1100 nm, the deviation of 60% might be also acceptable.

5.2 Comparison of the mean particle size of polystyrene latex beads from DLS and NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer.

The results of the mean particle size of polystyrene latex beads from DLS and NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer are listed in table 6.

Table 6. Mean particle size (Vol.) of polystyrene latex beads measured by DLS and NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer

Size of Polystyrene latex beads	PDI DLS		Mean particle size (D50, nm) from DLS (% deviation from standard)		Mean particle size (nm) NanoCuvette™ S and SpectroWorks™ together with UV-vis spectrophotometer (% deviation from standard)	
	Run 1	Run 2	Run 1	Run 2	Run 1	Run 2
100 nm	0.026	0.058	109 (9.0%)	107 (7.0%)	115.00 (15.0%)	148.90 (48.9%)
460 nm	0.025	0.048	442 (3.9%)	431 (6.3%)	414.10 (10.0%)	413.04 (10.2%)
600 nm	0.040	0.086	662 (10.3%)	633 (5.5%)	546.26 (9.0%)	549.85 (8.4%)
800 nm	0.070	0.102	704 (12.0%)	730 (8.8%)	633.44 (20.8%)	637.50 (20.3%)
1100nm	0.147	0.040	873 (20.6%)	953 (13.4%)	1097.66 (0.2%)	1083.85 (1.5%)
3000 nm	0.365	0.334	1690 (43.7%)	1480 (50.7%)	3637.24 (21.4%)	3050.95 (1.9%)

For the results from DLS measurements, the deviation from the standard increases with the increase of the particle size of the sample. This phenomenon is also in line with the increase of PDI with the increase in the particle size of the sample and might be due to more sedimentation of particles when



the particle size became larger. The results from replicate measurements done using NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer are consistent except for 100 nm (Table 6). Furthermore, particle size measurements done using NanoCuvette™ S with SpectroWorks™ gives the deviation less than 22% for all measurements except for replicate 100 nm measurements, whereas the deviation range from 3.9% to 50.7% for DLS measurements from the theoretical concentration. It should also be noted that the polystyrene latex beads provided by sigma has the standard deviation ranging from 5%-15%.

6 Conclusion

NanoCuvette™ S and SpectroWorks™ together with UV-Vis spectrophotometer is considered reliable to evaluate the concentration and particle size of the sample. The concentrations of all samples (except for polystyrene latex beads 1100 nm) from SpectroWorks™ have deviation less than 30% compared to the diluted samples from the stock provided by supplier. Furthermore, the particle size results obtained from SpectroWorks™ using measurements in UV-Vis spectrophotometer and NanoCuvette™ S have less than 22% deviation (except for run 2 of 100 nm sample) from the standard for all the size range.

MP CW



Change list

Issue	Change	Effective date
1	New document.	30-Jul-2021



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

**TESTING AND VALIDAITON OF NANOCUVETTE™ S, SPECTROWORKS™
TOGETHER WITH UV-VIS SPECTROPHOTOMETER**

Date of Issue: 30-Jul-2021

Supersedes: -

Appendix 1:

Results from NanoCuvette™ S and SpectroWork™ together with UV-Vis spectrophotometer (24 pages)

Summary saved

Experiment setup

Item ID:	f5d7d373
Box code:	NCONES
Cuvette no.:	01
Model:	NanoCuvette™ S
Created:	2021-07-26 11:11:42
Reference:	Water (n = 1.333215 nD)

Results


Refractive Index:	1.33612794 nD
Reference fit quality:	39.4927 %
Sample fit quality:	31.7062 %
Mean Particle Diameter:	115.00 nm
Particle Concentration:	0.0099 % (Vol.)

Sample attributes

	Sample name	Latex Beads 100nm 0.01%
	Protocol name	
	Analyte	
	Solvent	DI

Plots

Spectra



Size distribution

Here you can add notes about your sample.

Notes

Summary

Experiment setup

Item ID:	cb18820b
Box code:	NCONES
Cuvette no.:	01
Model:	NanoCuvette™ S
Created:	2021-07-26 11:26:48
Reference:	Water (n = 1.333215 nD)

Results


- Refractive Index: 1.33321537 nD
- Reference fit quality: 39.4927 %
- Sample fit quality: 25.9088 %
- Mean Particle Diameter: 148.90 nm
- Particle Concentration: 0.0081 % (Vol.)

Sample attributes

 Sample name	Latex beads 100nm 0.01% D1
 Protocol name	
 Analyte	
 Solvent	

Plots

Spectra 



Notes

Here you can add notes about your sample.

Summary

Experiment setup

Item ID:	f5d7d373
Box code:	NOONES
Cuvette no.:	01
Model:	NanoCuvette™ S
Created:	2021-07-26 11:11:42
Reference:	Water ($n = 1.333215 \text{ nD}$)


Results

Refractive Index:	1.33612794 nD
Reference fit quality:	39.4927 %
Sample fit quality:	31.7062 %
Mean Particle Diameter:	115.00 nm
Particle Concentration:	0.0099 % (Vol.)

Sample attributes

 Sample name	Larex Beads 100nm 0.01%
 Protocol name	
 Analyte	
 Solvent	DI

Plots



Notes

Here you can add notes about your sample.

Summary

Experiment setup

Item ID:	cb18820b
Box code:	NCONES
Cuvette no.:	01
Model:	NanoCuvette™ S
Created:	2021-07-26 11:26:48
Reference:	Water ($n = 1.333215 \text{ nD}$)


Results

- Refractive Index: 1.33321537 nD
- Reference fit quality: 39.4927 %
- Sample fit quality: 25.9088 %
- Mean Particle Diameter: 148.90 nm
- Particle Concentration: 0.0081 % (Vol.)

Sample attributes

 Sample name	Latex beads 100nm 0.01% Dl.
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	79e0fb8
Box code:	NCONES
Cuvette no.:	02
Model:	NanoCuvette™ S
Created:	2021-07-26 12:17:42
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: 1.36491180 nD
- Reference fit quality: 35.5789 %
- Sample fit quality: 18.5744 %
- Mean Particle Diameter: 275.73 nm
- Particle Concentration: 0.0038 % (Vol.)

Sample attributes

 Sample name	Latex Beads 460 nm 0.003%
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

not used discountinuity factor so concentration precise not size

Summary

Experiment setup

Item ID:	f23d4cb2
Box code:	NCONES
Cuvette no.:	02
Model:	NanoCuvette™ S
Created:	2021-07-26 12:45:42
Reference:	Water (n = 1.333215 nD)

Results


- Refractive Index: 1.36487618 nD
- Reference fit quality: 35.5789 %
- Sample fit quality: 17.1269 %
- Mean Particle Diameter: 284.26 nm
- Particle Concentration: 0.0038 % (Vol.)

Sample attributes

 Sample name	Latex Beads 460nm 0.003% D
 Protocol name	
 Analyte	
 Solvent	

Plots

Spectra 



Notes

Not used discontinuity factor so precise concentration

Summary saved

Experiment setup

Item ID: e6bd2474
 Box code: NCONES
 Cuvette no.: 02
 Model: NanoCuvette™ S
 Created: 2021-07-26 12:30:26
 Reference: Water ($n = 1.333215 \text{ nD}$)


Results

- Refractive Index: 1.36491180 nD
- Reference fit quality: 35.5789 %
- Sample fit quality: 18.6321 %
- Mean Particle Diameter: 414.10 nm
- Particle Concentration: 0.0082 % (Vol.)

Sample attributes

 Sample name	Latex Beads 460 nm 0.003%
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

Size precise Discontinuity used

Summary

Experiment setup

Item ID:	b6121334
Box code:	NCONES
Cuvette no.:	02
Model:	NanoCuvette™ S
Created:	2021-07-26 12:41:19
Reference:	Water (n = 1.3333215 nD)


Results

- Refractive Index: 1.36491179 nD
- Reference fit quality: 35.5789 %
- Sample fit quality: 17.1495 %
- Mean Particle Diameter: 413.04 nm
- Particle Concentration: 0.0085 % (Vol.)

Sample attributes

 Sample name	Latex Beads 460nm 0.003%
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

Size Precise
Discontinuity factor used

Summary saved

Experiment setup

Item ID: 2aaace06
Box code: NCONES
Cuvette no.: 03
Model: NanoCuvette™ S
Created: 2021-07-26 13:04:19
Reference: Water ($n = 1.333215 \text{ nD}$)

Results


Refractive Index:	1.35779775 nD
Reference fit quality:	54.7016 %
Sample fit quality:	28.0487 %
Mean Particle Diameter:	549.26 nm
Particle Concentration:	0.0026 % (Vol.)

Sample attributes

Sample name	Latex Beads 600nm 0.003%
Protocol name	
Analyte	
Solvent	

Plots

Spectra Size distribution



Notes

Here you can add notes about your sample.

Summary

Experiment setup

Item ID:	29668541
Box code:	NCONES
Cuvette no.:	03
Model:	NanoCuvette™ S
Created:	2021-07-26 13:08:53
Reference:	Water ($n = 1.333215$ nD)


Results

- Refractive Index: 1.35936624 nD
- Reference fit quality: 54.7016 %
- Sample fit quality: 28.8134 %
- Mean Particle Diameter: 549.85 nm
- Particle Concentration: 0.0026 % (Vol.)

Sample attributes

 Sample name	Latex Beads 600nm 0.003% D
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	2aaece06
Box code:	NCONES
Cuvette no.:	03
Model:	NanoCuvette™ S
Created:	2021-07-26 13:04:19
Reference:	Water ($n = 1.333215 \text{ nD}$)


Results

- Refractive Index: 1.35779775 nD
- Reference fit quality: 54.7016 %
- Sample fit quality: 28.0487 %
- Mean Particle Diameter: 549.26 nm
- Particle Concentration: 0.0026 % (Vol.)

Sample attributes

 Sample name	Latex Beads 600nm 0.003%
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	29f68541
Box code:	NCONES
Cuvette no.:	03
Model:	NanoCuvette™ S
Created:	2021-07-26 13:08:53
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: 1.35936624 nD
- Reference fit quality: 54.7016 %
- Sample fit quality: 28.8134 %
- Mean Particle Diameter: 549.85 nm
- Particle Concentration: 0.0026 % (Vol.)

Sample attributes

 Sample name	Latex Beads 600nm 0.003% D
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	1d43d47f
Box code:	NCONES
Cuvette no.:	04
Model:	NanoCuvette™ S
Created:	2021-07-26 13:33:07
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: 1.34305982 nD
- Reference fit quality: 55.9293 %
- Sample fit quality: 24.9135 %
- Mean Particle Diameter: 633.44 nm
- Particle Concentration: 0.0030 % (Vol.)

Sample attributes

<input type="checkbox"/> Sample name	Latex Beads 800nm 0.003%
<input type="checkbox"/> Protocol name	
<input type="checkbox"/> Analyte	
<input type="checkbox"/> Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary

Experiment setup

Item ID:	2b515f9e
Box code:	NCONES
Cuvette no.:	04
Model:	NanoCuvette™ S
Created:	2021-07-26 13:36:32
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: 1.34309282 nD
- Reference fit quality: 55.9293 %
- Sample fit quality: 24.7994 %
- Mean Particle Diameter: 637.50 nm
- Particle Concentration: 0.0031 % (Vol.)

Sample attributes

<input type="checkbox"/> Sample name	Latex Beads 800nm 0.003%
<input type="checkbox"/> Protocol name	
<input type="checkbox"/> Analyte	
<input type="checkbox"/> Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary

Experiment setup

Item ID:	1d43d47f
Box code:	NCONES
Cuvette no.:	04
Model:	NanoCuvette™ S
Created:	2021-07-26 13:33:07
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: 1.34305982 nD
- Reference fit quality: 55.9293 %
- Sample fit quality: 24.9135 %
- Mean Particle Diameter: 633.44 nm
- Particle Concentration: 0.0030 % (Vol.)

Sample attributes

<input checked="" type="checkbox"/> Sample name	Latex Beads 800nm 0.003%
<input checked="" type="checkbox"/> Protocol name	
<input checked="" type="checkbox"/> Analyte	
<input checked="" type="checkbox"/> Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	2b515f9e
Box code:	NCONES
Cuvette no.:	04
Model:	NanoCuvette™ S
Created:	2021-07-26 13:36:32
Reference:	Water (n = 1.333215 nD)


Results

Refractive Index:	1.34309282 nD
Reference fit quality:	55.9293 %
Sample fit quality:	24.7994 %
Mean Particle Diameter:	637.50 nm
Particle Concentration:	0.0031 % (Vol.)

Sample attributes

 Sample name	Latex Beads 800nm 0.003% D
 Protocol name	
 Analyte	
 Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	3ece3e8d
Box code:	NCONES
Cuvette no.:	05
Model:	NanoCuvette™ S
Created:	2021-07-26 13:54:00
Reference:	Water (n = 1.333215 nD)


Results

Refractive Index:	N/A
Reference fit quality:	N/A
Sample fit quality:	N/A
Mean Particle Diameter:	1097.66 nm
Particle Concentration:	0.0016 % (Vol.)

Sample attributes

<input type="checkbox"/> Sample name	Latex Beads 1100nm 0.001%
<input type="checkbox"/> Protocol name	
<input type="checkbox"/> Analyte	
<input type="checkbox"/> Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	0f5caa60
Box code:	NCONES
Cuvette no.:	05
Model:	NanoCuvette™ S
Created:	2021-07-26 14:00:18
Reference:	Water (n = 1.333215 nD)


Results

- Refractive index: N/A
- Reference fit quality: N/A
- Sample fit quality: N/A
- Mean Particle Diameter: 1083.85 nm
- Particle Concentration: 0.0016 % (Vol.)

Sample attributes

<input type="checkbox"/> Sample name	Latex Beads 110nm 0.001%
<input type="checkbox"/> Protocol name	
<input type="checkbox"/> Analyte	
<input type="checkbox"/> Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	3ece3e8d
Box code:	NCONES
Cuvette no.:	05
Model:	NanoCuvette™ S
Created:	2021-07-26 13:54:00
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: N/A
- Reference fit quality: N/A
- Sample fit quality: N/A
- Mean Particle Diameter: 1097.66 nm
- Particle Concentration: 0.0016 % (Vol.)

Sample attributes

<input checked="" type="checkbox"/>	Sample name	Latex Beads 1100nm 0.001%
<input type="checkbox"/>	Protocol name	
<input type="checkbox"/>	Analyte	
<input type="checkbox"/>	Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	0f5caa60
Box code:	NCONES
Cuvette no.:	05
Model:	NanoCuvette™ S
Created:	2021-07-26 14:00:18
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: N/A
- Reference fit quality: N/A
- Sample fit quality: N/A
- Mean Particle Diameter: 1083.85 nm
- Particle Concentration: 0.0016 % (Vol.)

Sample attributes

<input type="checkbox"/> Sample name	Latex Beads 1100nm 0.001%
<input type="checkbox"/> Protocol name	
<input type="checkbox"/> Analyte	
<input type="checkbox"/> Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	805ae0d2
Box code:	NCONES
Cuvette no.:	06
Model:	NanoCuvette™ S
Created:	2021-07-26 14:18:16
Reference:	Water (n = 1.333215 nD)

Results


Refractive Index:	1.88416545 nD
Reference fit quality:	35.7368 %
Sample fit quality:	23.5322 %
Mean Particle Diameter:	3637.24 nm
Particle Concentration:	0.0008 % (Vol.)

Sample attributes

▷ Sample name	Latex Beads 3000nm 0.001%
▷ Protocol name	
▷ Analyte	
▷ Solvent	

Plots

Spectra



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	551fa119
Box code:	NCONES
Cuvette no.:	06
Model:	NanoCuvette™ S
Created:	2021-07-26 14:22:39
Reference:	Water ($n = 1.333215 \text{ nD}$)

Results


- Refractive Index: 1.57637491 nD
- Reference fit quality: 35.7368 %
- Sample fit quality: 59.9385 %
- Mean Particle Diameter: 3050.95 nm
- Particle Concentration: 0.0007 % (Vol.)

Sample attributes

 Sample name	Latex Beads 3000nm 0.001%
 Protocol name	
 Analyte	
 Solvent	

Plots

Spectra  Size distribution



Notes

Here you can add notes about your sample.

Summary

23

Experiment setup

Item ID:	805ae0d2
Box code:	NCONES
Cuvette no.:	06
Model:	NanoCuvette™ S
Created:	2021-07-26 14:18:16
Reference:	Water (n = 1.333215 nD)


Results

Refractive Index:	1.88416545 nD
Reference fit quality:	35.7368 %
Sample fit quality:	23.5322 %
Mean Particle Diameter:	3637.24 nm
Particle Concentration:	0.0008 % (Vol.)

Sample attributes

Sample name	Latex Beads 3000nm 0.001%
Protocol name	
Analyte	
Solvent	

Plots



Notes

Here you can add notes about your sample.

Summary saved

Experiment setup

Item ID:	551fa119
Box code:	NCONES
Cuvette no.:	06
Model:	NanoCuvette™ S
Created:	2021-07-26 14:22:39
Reference:	Water (n = 1.333215 nD)


Results

- Refractive Index: 1.57637491 nD
- Reference fit quality: 35.7368 %
- Sample fit quality: 59.9385 %
- Mean Particle Diameter: 3050.95 nm
- Particle Concentration: 0.0007 % (Vol.)

Sample attributes

<input type="checkbox"/> Sample name	Latex Beads 3000nm 0.001%
<input type="checkbox"/> Protocol name	
<input type="checkbox"/> Analyte	
<input type="checkbox"/> Solvent	

Plots



Notes

Here you can add notes about your sample.



PARTICLE
ANALYTICAL

EXPERTS IN SIZE AND CRYSTALS

**TESTING AND VALIDAITON OF NANOCUVETTE™ S, SPECTROWORKS™
TOGETHER WITH UV-VIS SPECTROPHOTOMETER**

Date of Issue: 30-Jul-2021

Supersedes: -

Appendix 2:

Results from dynamic light scattering (60 pages)



CERTIFICATE OF ANALYSIS

Customer: CphNano
Material tested: Polystyrene Latex beads 100nm
Batch: 2021/07/23
Internal number: 32027

Analytical technique: Dynamic Light Scattering
Method of analysis: Analyseplan
Internal quality level: GMP

	D _{10%} (nm)**	D _{50%} (nm)**	D _{90%} (nm)**	Z-average*	PDI*
Run 1	79.4	109	159	118	0.026
Run 2	74.7	107	163	118.1	0.058
Average	76.9	108	161	118.1	0.042
Specification	-	-	-	-	-
Evaluation	-	-	-	-	-

*The parameter is based on the intensity size distribution

**The parameter is based on the volume size distribution

Written by: Wenbo Wang

Date: 28 JULI 2021

Reviewed by: Wenbo Wang 28 JULI 2021^①

Date: 28 JULI 2021

Approved by: Wenbo Wang 28 JULI 2021^①
QC

Date: 28 JULI 2021

The validity of the method is the responsibility of the sponsor
Quality agreement not in place

① wrong stamp up 28-Jul-2021

Size Distribution Report by Volume


Sample Name: 32027 batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 08:48:19

Userid: brj

D10%(V): 78.9 nm			D50%(V): 109 nm	D90%(V): 160 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and Pdl are based on the intensity distribution
Peak 1:	114.3	100.0	30.95			
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	259.1	



General Notes: 32027 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator: 27 JULI 2021 BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved:

28 JULI 2021 WJP

File name: 32027
Record Number: 2
27 jul 2021 14:08




Sample Name: 32027 batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 08:50:33

Userid: brj

D10%(V): 79.8 nm			D50%(V): 109 nm	D90%(V): 158 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution
Peak 1:	114.1	100.0	29.46	118.1	0.014	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	260.6	



General Notes: 32027 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator:
 Malvern Panalytical
 www.malvernpanalytical.com

27 JULI 2021
 BRJ

Zetasizer Ver. 7.13
 Serial Number : MAL500686

Approved:

28 JULI 2021

WJ

File name: 32027
 Record Number: C
 27 jul 2021 14:08

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32027 batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 08:52:46

Userid: brj

D10%(V): 79.5 nm			D50%(V): 109 nm	D90%(V): 159 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution
Peak 1:	114.0	100.0	29.84			
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	258.6	



General Notes: 32027 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	
		Size Runs:	13

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021

Approved:

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32027
Record Number: 4
27 jul 2021 14:08




Sample Name: 32027 batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 08:59:41

Userid: brj

D10%(V): 74.3 nm			D50%(V): 107 nm	D90%(V): 164 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and Pdl are based on the intensity distribution
Peak 1:	113.8	100.0	34.99			
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	242.7	



General Notes: 32027 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	14

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021
BRJ

Approved:

28 JULI 2021
WJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32027
Record Number: f
27 jul 2021 14:08



Sample Name: 32027 batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 09:02:05

Userid: brj

D10%(V): 72.4 nm			D50%(V): 105 nm	D90%(V): 164 nm
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):
Peak 1:	112.3	100.0	35.93	PDI: 0.068
Peak 2:	0.000	0.0	0.000	
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 243.9

Note: Z-average and PDI are based on the intensity distribution



General Notes: 32027 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10 Extend duration for large particles: False
		Size Runs:	14

Operator:
 Malvern Panalytical
 www.malvernpanalytical.com

27 JULI 2021

Approved:

Zetasizer Ver. 7.13
 Serial Number : MAL500686

28 JULI 2021

WFO

File name: 32027
 Record Number: f
 27 jul 2021 14:09

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32027 batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 09:04:29

Userid: brj

	D10%(V): 78.3 nm	D50%(V): 109 nm	D90%(V): 162 nm
Peak 1:	Diam. (nm) 115.1	% Volume 100.0	Width (nm) 32.42
Peak 2:	0.000	0.0	Z-Average (d.nm): 118.7 PDI: 0.048
Peak 3:	0.000	0.0	Count Rate (kcps): 243.2



General Notes: 32027 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	14

Operator:

27 JULI 2021 BRJ

Approved:

28 JULI 2021

BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32027
Record Number: 7
27 jul 2021 14:09



Size Distribution Report by Volume


Sample Name: 32027 batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 09:39:35

Userid: brj

D10%(V): 79.4 nm			D50%(V): 109 nm	D90%(V): 159 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	118.0	Note: Z-average and Pdl are based on the intensity distribution
Peak 1:	114.1	100.0	30.09	PdI:	0.026	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	259.1	



General Notes: Average result created from record number(s): 2 3 4

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator:

27 JULI 2021 *[Signature]*

Approved:

28 JULI 2021 *[Signature]*

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32027
Record Number: 8
27 jul 2021 14:09:



CERTIFICATE OF ANALYSIS

Customer: CphNano
Material tested: Polystyrene Latex beads 460nm
Batch: 2021/07/23
Internal number: 32028

Analytical technique: Dynamic Light Scattering
Method of analysis: Analyseplan
Internal quality level: GMP

	D _{10%} (nm)**	D _{50%} (nm)**	D _{90%} (nm)**	Z-average*	PDI*
Run 1	320	442	605	422.2	0.025
Run 2	316	431	587	413.6	0.048
Average	318	436	597	417.9	0.037
Specification	-	-	-	-	-
Evaluation	-	-	-	-	-

*The parameter is based on the intensity size distribution

**The parameter is based on the volume size distribution

Written by: Wenbo Wang

Date: 28 JULI 2021

Reviewed by: Wenbo Wang

Date: 28 JULI 2021

Approved by: Wenbo Wang
QC

Date: 28 JULI 2021

The validity of the method is the responsibility of the sponsor
Quality agreement not in place

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32028 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27.juli 2021 10:03:54

Userid: brj

	D10%(V): 313 nm	D50%(V): 428 nm	D90%(V): 585 nm
Peak 1:	Diam. (nm) 437.0	% Volume 100.0	Width (nm) 96.69
Peak 2:	0.000	0.0	Z-Average (d.nm): 411.4 PDI: 0.035
Peak 3:	0.000	0.0	Count Rate (kcps): 423.7

Note: Z-average and PDI are based on the intensity distribution



General Notes: 32028 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	12

Operator: 27 JULI 2021 BRJ
Malvern Panalytical
www.malvernpanalytical.com

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved: WJ

28 JULI 2021

File name: 32028
Record Number: 1
27 jul 2021 14:11

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32028 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:05:57

Userid: brj

D10%(V): 328	nm	D50%(V): 448	nm	D90%(V): 604	nm
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	430.2
Peak 1:	457.1	100.0	99.23	PDI:	0.036
Peak 2:	0.000	0.0	0.000	Note: Z-average and PDI are based on the intensity distribution	
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 427.0	



General Notes: 32028 Run 1

Cell Description: Disposable sizing cuvette

Material RI: 1.59

Dispersant Name: Water

Material Absorbtion: 0.010

Dispersant RI: 1.330

Analysis Model: General Purpose

Viscosity (cP): 0.8872

Lower Size Threshold: 0.050

Temperature (°C): 25.0

Upper Size Threshold: 0.010

Equilibration Time Set (min): 180

Size range: 0.6000 to 6000 nm

Number of measurements: 3

Size Measure Delay (s): 0

Auto Position Enabled: False

Measurement Position (mm): 4.65

Auto Attenuate Enabled: True

Attenuator: 7

Auto Size Measurement Time: True

Duration (s): 10 Extend duration for large particles: False

Size Runs: 12

Operator: 27 JULI 2021 BRJ

Approved: *WJ*

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32028
Record Number: 2
27 jul 2021 14:11




Sample Name: 32028 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:08:00

Userid: brj

D10%(V): 321	nm	D50%(V): 451	nm	D90%(V): 630	nm
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	425.0
Peak 1:	463.8	100.0	113.0	PdI:	0.004
Peak 2:	0.000	0.0	0.000		
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	428.7



General Notes: 32028 Run 1

Cell Description: Disposable sizing cuvette

Material RI: 1.59

Dispersant Name: Water

Material Absorbtion: 0.010

Dispersant RI: 1.330

Analysis Model: General Purpose

Viscosity (cP): 0.8872

Lower Size Threshold: 0.050

Temperature (°C): 25.0

Upper Size Threshold: 0.010

Equilibration Time Set (min): 180

Size range: 0.6000 to 6000 nm

Number of measurements: 3

Size Measure Delay (s): 0

Auto Position Enabled: False

Measurement Position (mm): 4.65

Auto Attenuate Enabled: True

Attenuator: 7

Auto Size Measurement Time: True

Duration (s): 10 Extend duration for large particles: False

Size Runs: 12

Operator: brjmalvern Panalytical
www.malvernpanalytical.comApproved: WJZetasizer Ver. 7.13
Serial Number : MAL500686File name: 32028
Record Number: 3
27 jul 2021 14:11

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32028 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:15:25

Userid: brj

D10%(V): 319 nm			D50%(V): 436 nm			D90%(V): 592 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	420.0	Note: Z-average and PDI are based on the intensity distribution		
Peak 1:	445.0	100.0	97.11	PDI:	0.047			
Peak 2:	0.000	0.0	0.000					
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 420.3				



General Notes: 32028 Run 2

Cell Description: Disposable sizing cuvette

Material RI: 1.59

Dispersant Name: Water

Material Absorbtion: 0.010

Dispersant RI: 1.330

Analysis Model: General Purpose

Viscosity (cP): 0.8872

Lower Size Threshold: 0.050

Temperature (°C): 25.0

Upper Size Threshold: 0.010

Equilibration Time Set (min): 180

Size range: 0.6000 to 6000 nm

Number of measurements: 3

Size Measure Delay (s): 0

Auto Position Enabled: False

Measurement Position (mm): 4.65

Auto Attenuate Enabled

True

Attenuator: 7

Auto Size Measurement Time:

True

Duration (s): 10

Extend duration for large particles:

False

Size Runs: 12

Operator:

27 JULI 2021 BRJ

Approved:

28 JULI 2021

WPA

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32028
Record Number: 4
27 jul 2021 14:11

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32028 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:17:28

Userid: brj

	D10%(V): 310 nm	D50%(V): 419 nm	D90%(V): 570 nm
Peak 1:	Diam. (nm) 428.0	% Volume 100.0	Width (nm) 92.11
Peak 2:	0.000	0.0	Z-Average (d.nm): 397.4 PDI: 0.036
Peak 3:	0.000	0.0	Width (nm) 0.000 Count Rate (kcps): 418.5

Note: Z-average and PDI are based on the intensity distribution



General Notes: 32028 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	12

Operator: 27 JULI 2021 BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved:

28 JULI 2021

File name: 32028
Record Number: 5
27 jul 2021 14:11

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32028 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:19:31

Userid: brj

D10%(V):	320 nm	D50%(V):	438 nm	D90%(V):	596 nm
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	423.5
Peak 1:	447.5	100.0	98.21	PdI:	0.061
Peak 2:	0.000	0.0	0.000		
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	419.4

Size Distribution by Volume



Record 6: 32028 Batch 2021/07/23

General Notes: 32028 Run 2

Cell Description: Disposable sizing cuvette

Material RI: 1.59

Dispersant Name: Water

Material Absorbtion: 0.010

Dispersant RI: 1.330

Analysis Model: General Purpose

Viscosity (cP): 0.8872

Lower Size Threshold: 0.050

Temperature (°C): 25.0

Upper Size Threshold: 0.010

Equilibration Time Set (min): 180

Size range: 0.6000 to 6000 nm

Number of measurements: 3

Size Measure Delay (s): 0

Auto Position Enabled: False

Measurement Position (mm): 4.65

Auto Attenuate Enabled: True

Attenuator: 7

Auto Size Measurement Time:

True

Duration (s): 10

Extend duration for large particles: False

Size Runs: 12

Operator: *BRJ*

27 JULI 2021 BRJ

Approved:

28 JULI 2021

WPA

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32028
Record Number: f
27 jul 2021 14:11

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32028 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:20:04

Userid: brj

D10%(V): 320 nm			D50%(V): 442 nm	D90%(V): 605 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	422.2
Peak 1:	452.6	100.0	103.8	PDI:	0.025
Peak 2:	0.000	0.0	0.000		
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	423.7

Note: Z-average and Pdl are based on the intensity distribution



General Notes: Average result created from record number(s): 1 2 3

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	12

Operator:

27 JULI 2021 BRJ

Approved:

28 JULI 2021

W.M.

File name: 32028
Record Number: 7
27 jul 2021 14:11

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32028 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:20:16

Userid: brj

	D10%(V): 316 nm	D50%(V): 431 nm	D90%(V): 587 nm
	Diam. (nm)	% Volume	Width (nm)
Peak 1:	440.2	100.0	96.24
Peak 2:	0.000	0.0	0.000
Peak 3:	0.000	0.0	0.000
		Z-Average (d.nm):	413.6
		PdI:	0.048
		Count Rate (kcps):	420.3

Note: Z-average and PdI are based on the intensity distribution



General Notes: Average result created from record number(s): 4 5 6

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	12

Operator:
alvern Panalytical
www.alvernpanalytical.com

27 JULI 2021

Approved:

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32028
Record Number: 8
27 jul 2021 14:11

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32028 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27.juli 2021 10:20:28

Userid: brj

D10%(V):	318	nm	D50%(V):	436	nm	D90%(V):	597	nm
		Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	417.9		
Peak 1:	446.4		100.0	100.3	PdI:	0.037	<i>Note: Z-average and PdI are based on the intensity distribution</i>	
Peak 2:	0.000		0.0	0.000				
Peak 3:	0.000		0.0	0.000	Count Rate (kcps):	423.7		



General Notes: Average result created from record number(s): 1 2 3 4 5 6

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000	nm
			Number of measurements:	3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	7
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	12

Operator: 27 JULI 2021

malvern Panalytical
www.malvernpanalytical.com

Approved:

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32028
Record Number: 9
27 jul 2021 14:11



CERTIFICATE OF ANALYSIS

Customer: CphNano
Material tested: Polystyrene Latex beads 600nm
Batch: 2021/07/23
Internal number: 32029

Analytical technique: Dynamic Light Scattering
Method of analysis: Analyseplan
Internal quality level: GMP

	D _{10%} (nm)**	D _{50%} (nm)**	D _{90%} (nm)**	Z-average*	PDI*
Run 1	449	662	948	569.5	0.04
Run 2	470	633	831	596.4	0.086
Average	462	646	899	582.9	0.063
Specification	-	-	-	-	-
Evaluation	-	-	-	-	-

*The parameter is based on the intensity size distribution

**The parameter is based on the volume size distribution

Written by: Bente Ansga Jøn

Date: 28 JULI 2021

Reviewed by: Wenbo Wang

Date: 28 JULI 2021

Approved by: Wenbo Wang
QC

Date: 28 JULI 2021

The validity of the method is the responsibility of the sponsor
Quality agreement not in place




Sample Name: 32029 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:32:39

Userid: brj

D10%(V): 414 nm			D50%(V): 705 nm			D90%(V): 1090 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	538.3	Note: Z-average and PDI are based on the intensity distribution	
Peak 1:	729.3	100.0	250.1	PDI:	0.090		
Peak 2:	0.000	0.0	0.000				
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		188.5	



General Notes: 32029 Run 1

Cell Description: Disposable sizing cuvette

Material RI: 1.59

Dispersant Name: Water

Material Absorbtion: 0.010

Dispersant RI: 1.330

Analysis Model: General Purpose

Viscosity (cP): 0.8872

Lower Size Threshold: 0.050

Temperature (°C): 25.0

Upper Size Threshold: 0.010

Equilibration Time Set (min): 180

Size range: 0.6000 to 6000 nm

Number of measurements: 3

Size Measure Delay (s): 0

Auto Position Enabled: False

Measurement Position (mm): 4.65

Auto Attenuate Enabled: True

Attenuator: 7

Auto Size Measurement Time: True

Duration (s): 10 Extend duration for large particles: False

Size Runs: 15

Operator: 27 JULI 2021 BRJ
 Malvern Panalytical
 www.malvernpanalytical.com

 Zetasizer Ver. 7.13
 Serial Number : MAL500686
Approved: 28 JULI 2021 WLN
 File name: 32027
 Record Number:
 27 jul 2021 14:10




Sample Name: 32029 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:35:13

Userid: brj

D10%(V): 455	nm	D50%(V): 637	nm	D90%(V): 873	nm
Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and Pdl are based on the intensity distribution
Peak 1: 647.4	100.0	149.7	573.1	0.024	
Peak 2: 0.000	0.0	0.000			
Peak 3: 0.000	0.0	0.000	Count Rate (kcps):	188.1	



General Notes: 32029 Run 1

Cell Description:	Disposable sizing cuvette				
Material RI:	1.59	Dispersant Name:	Water		
Material Absorbtion:	0.010	Dispersant RI:	1.330		
Analysis Model:	General Purpose	Viscosity (cP):	0.8872		
Lower Size Threshold:	0.050	Temperature (°C):	25.0		
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180		
Size range:	0.6000 to 6000 nm	Number of measurements:	3		
Auto Position Enabled:	False	Size Measure Delay (s):	0		
Auto Attenuate Enabled	True	Measurement Position (mm):	4.65		
Auto Size Measurement Time:	True	Attenuator:	7	Extend duration for large particles:	False
		Duration (s):	10	Size Runs:	15

Operator: 27 JULI 2021 BRJ
 alvern Panalytical
 www.alvernpanalytical.com

Zetasizer Ver. 7.13
 Serial Number : MAL500686

Approved:

28 JULI 2021 WJ

File name: 32027
 Record Number: 27 jul 2021 14:10:

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32029 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:37:46

Userid: brj

D10%(V): 474 nm			D50%(V): 661 nm			D90%(V): 902 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution	
Peak 1:	671.4	100.0	151.9			597.1	
Peak 2:	0.000	0.0	0.000			0.005	
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		188.9	



General Notes: 32029 Run 1

Cell Description: Disposable sizing cuvette

Material RI: 1.59

Dispersant Name: Water

Material Absorbtion: 0.010

Dispersant RI: 1.330

Analysis Model: General Purpose

Viscosity (cP): 0.8872

Lower Size Threshold: 0.050

Temperature (°C): 25.0

Upper Size Threshold: 0.010

Equilibration Time Set (min): 180

Size range: 0.6000 to 6000 nm

Number of measurements: 3

Size Measure Delay (s): 0

Auto Position Enabled: False

Measurement Position (mm): 4.65

Auto Attenuate Enabled: True

Attenuator: 7

Auto Size Measurement Time:

True

Duration (s): 10

Extend duration for large particles:

False

Size Runs: 15

Operator: alvern Panalytical

27 JULI 2021 BRJ

Approved:

28 JULI 2021

WJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32029
Record Number: 13
27 jul 2021 14:10



Sample Name: 32029 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:44:42

Userid: brj

	D10%(V): 475 nm	D50%(V): 636 nm	D90%(V): 829 nm
Peak 1:	Diam. (nm) 644.1	% Volume 100.0	Width (nm) 127.2
Peak 2:	0.000	0.0	Z-Average (d.nm): 607.3 PDI: 0.134
Peak 3:	0.000	0.0	Count Rate (kcps): 175.4

Note: Z-average and PDI are based on the intensity distribution



General Notes: 32029 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	15

Operator: 27 JULI 2021 BRJ
 malvern Panalytical
 www.malvernpanalytical.com

Zetasizer Ver. 7.13
 Serial Number : MAL500686

Approved: 28 JULI 2021 WJ
 File name: 32027
 Record Number:
 27 jul 2021 14:10:




Sample Name: 32029 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:47:15

Userid: brj

D10%(V): 474 nm			D50%(V): 633 nm			D90%(V): 824 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PdI:	Note: Z-average and PdI are based on the intensity distribution		
Peak 1:	641.4	100.0	126.4			600.4		
Peak 2:	0.000	0.0	0.000					
Peak 3:	0.000	0.0	0.000			175.8		



General Notes: 32029 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	15

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021 BRJ

Approved:

Zetasizer Ver. 7.13
Serial Number : MAL500686

28 JULI 2021

WJ
File name: 32027
Record Number:
27 jul 2021 14:10

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32029 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:49:49

Userid: brj

	D10%(V): 460 nm	D50%(V): 631 nm	D90%(V): 839 nm
	Diam. (nm)	% Volume	Width (nm)
Peak 1:	638.3	100.0	136.4
Peak 2:	0.000	0.0	0.000
Peak 3:	0.000	0.0	0.000
		Z-Average (d.nm):	581.5
		PDI:	0.027
		Count Rate (kcps):	173.2



Record 16: 32029 Batch 2021/07/23

General Notes: 32029 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	15

Operator: 27 JULI 2021

Approved:

28 JULI 2021

WJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32029
Record Number:
27 jul 2021 14:10

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32029 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:51:35

Userid: brj

D10%(V): 449 nm			D50%(V): 662 nm	D90%(V): 948 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution
Peak 1:	682.7	100.0	192.8			
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	188.5	



General Notes: Average result created from record number(s): 11 12 13

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59	Dispersant Name:	Water	
Material Absorbtion:	0.010	Dispersant RI:	1.330	
Analysis Model:	General Purpose	Viscosity (cP):	0.8872	
Lower Size Threshold:	0.050	Temperature (°C):	25.0	
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180	
Size range:	0.6000 to 6000 nm	Number of measurements:	3	
		Size Measure Delay (s):	0	
Auto Position Enabled:	False	Measurement Position (mm):	4.65	
Auto Attenuate Enabled	True	Attenuator:	7	
Auto Size Measurement Time:	True	Duration (s):	10	Extend duration for large particles:
		Size Runs:	15	False

Operator:
Malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021
brj

Approved:

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32027
Record Number: 27 jul 2021 14:10:

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32029 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:51:46

Userid: brj

	D10%(V): 470 nm	D50%(V): 633 nm	D90%(V): 831 nm
Peak 1:	Diam. (nm) 641.3	% Volume 100.0	Width (nm) 130.1
Peak 2:	0.000	0.0	Z-Average (d.nm): 596.4 PDI: 0.086
Peak 3:	0.000	0.0	Count Rate (kcps): 175.4



General Notes: Average result created from record number(s): 14 15 16

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10 Extend duration for large particles: False
		Size Runs:	15

28 JULI 2021

Approved:

File name: 32027
Record Number:
27 jul 2021 14:10

Operator: 27 JULI 2021 BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32029 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 10:52:50

Userid: brj

D10%(V): 462 nm			D50%(V): 646 nm	D90%(V): 899 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	582.9	Note: Z-average and Pdl are based on the intensity distribution
Peak 1:	662.0	100.0	165.8	PdI:	0.063	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	188.5	

Size Distribution by Volume



Record 19: 32029 Batch 2021/07/23

General Notes: Average result created from record number(s): 11 12 13 14 15 16

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	7
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	15

Operator:

27 JULI 2021 BRJ

Approved:

28 JULI 2021

BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32027
Record Number: 27 jul 2021 14:10



CERTIFICATE OF ANALYSIS

Customer: CphNano
Material tested: Polystyrene Latex beads 800nm
Batch: 2021/07/23
Internal number: 32030

Analytical technique: Dynamic Light Scattering
Method of analysis: Analyseplan
Internal quality level: GMP

	D _{10%} (nm)**	D _{50%} (nm)**	D _{90%} (nm)**	Z-average*	PDI*
Run 1	496	704	948	621.9	0.07
Run 2	552	730	944	684.9	0.102
Average	529	718	946	653.4	0.086
Specification	-	-	-	-	-
Evaluation	-	-	-	-	-

*The parameter is based on the intensity size distribution

**The parameter is based on the volume size distribution

Written by: Boštjan Štefanec

Date: 28 JULI 2021

Reviewed by: Wenbo Wang

Date: 28 JULI 2021

Approved by: Wenbo Wang
QC

Date: 28 JULI 2021

The validity of the method is the responsibility of the sponsor
Quality agreement not in place

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:04:35

Userid: brj

D10%(V): 436 nm			D50%(V): 670 nm			D90%(V): 968 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution		
Peak 1:	686.9	100.0	196.5					
Peak 2:	0.000	0.0	0.000					
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 290.7				



General Notes: 32030 Run1

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000 nm	
			Number of measurements:	3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	8
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	13

Operator: 27 JULI 2021 BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved: 28 JULI 2021

WJ

File name: 32030
Record Number: 27 jul 2021 14:12




Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:06:48

Userid: brj

	D10%(V): 515 nm	Diam. (nm)	% Volume	Width (nm)	D50%(V): 700 nm	Z-Average (d.nm):	643.0	D90%(V): 933 nm
Peak 1:	711.2	100.0	149.1			PdI:	0.038	
Peak 2:	0.000	0.0	0.000					
Peak 3:	0.000	0.0	0.000		Count Rate (kcps):		285.5	



General Notes: 32030 Run1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	8
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator: 27 JULI 2021 BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved: 28 JULI 2021

WJ

File name: 32030
Record Number: 2
27 jul 2021 14:12:

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:09:01

Userid: brj

D10%(V): 553 nm			D50%(V): 735 nm	D90%(V): 951 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	675.1	Note: Z-average and PDI are based on the intensity distribution
Peak 1:	744.2	100.0	142.4	PdI:	0.008	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 288.2		



General Notes: 32030 Run1

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000 nm	Number of measurements: 3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	8
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	13

Operator: 27 JULI 2021

Approved: 28 JULI 2021

VWA

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32030
Record Number: C
27 jul 2021 14:12

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27.juli 2021 11:16:03

Userid: brj

D10%(V): 543 nm			D50%(V): 710 nm			D90%(V): 925 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	669.6	Note: Z-average and Pdl are based on the intensity distribution	
Peak 1:	720.1	100.0	134.0	PdI:	0.089		
Peak 2:	0.000	0.0	0.000				
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 275.1			



General Notes: 32030 Run 2

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000	nm
			Number of measurements:	3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	8
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	13

Operator:

27 JULI 2021 BRJ

Approved:

WJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32030
Record Number: 4
27 jul 2021 14:12

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27.juli 2021 11:18:16

Userid: brj

D10%(V): 565 nm			D50%(V): 748 nm			D90%(V): 958 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	697.1	Note: Z-average and PDI are based on the intensity distribution	
Peak 1:	755.4	100.0	139.9	PDI:	0.097		
Peak 2:	0.000	0.0	0.000				
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		274.4	



General Notes: 32030 Run 2

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000 nm	
			Number of measurements:	3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	8
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	13

Operator:

27 JULI 2021 BRJ

Approved:

28 JULI 2021 WYD

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32030
Record Number: 5
27 jul 2021 14:12

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:20:29

Userid: brj

D10%(V): 552 nm			D50%(V): 731 nm	D90%(V): 946 nm
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):
Peak 1:	739.8	100.0	140.2	PDI: 0.118
Peak 2:	0.000	0.0	0.000	
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 274.4



General Notes: 32030 Run 2

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000 nm	
			Number of measurements:	3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	8
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	13

Operator:

27 JULI 2021 BRJ

Approved:

28 JULI 2021 WY

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32030
Record Number: f
27 jul 2021 14:12:

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:21:24

Userid: brj

D10%(V): 496 nm			D50%(V): 704 nm	D90%(V): 948 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution
Peak 1:	714.1	100.0	166.1	621.9	0.070	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	290.7	



General Notes: Average result created from record number(s): 1 2 3

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	8
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator: *BRJ*

Malvern Panalytical
www.malvernpanalytical.com

Approved: *WY*

Zetasizer Ver. 7.13
Serial Number : MAL500686

28 JULI 2021

WY

File name: 32030
Record Number: 7
27 jul 2021 14:12

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32030 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:21:35

Userid: brj

D10%(V): 552 nm			D50%(V): 730 nm	D90%(V): 944 nm
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):
Peak 1:	738.4	100.0	138.8	684.9
Peak 2:	0.000	0.0	0.000	PdI: 0.102
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 275.1



General Notes: Average result created from record number(s): 4 5 6

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	8
Auto Size Measurement Time:	True	Duration (s):	10 Extend duration for large particles: False
		Size Runs:	13

Operator: 27 JULI 2021 BRJ
alvern Panalytical
www.alvernpanalytical.com

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved: 28 JULI 2021 WWD

File name: 32030
Record Number: 8
27 jul 2021 14:12

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32030 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:21:47

Userid: brj

	D10%(V): 529 nm	D50%(V): 718 nm	D90%(V): 946 nm
Peak 1:	Diam. (nm) 726.3	% Volume 100.0	Width (nm) 153.6
Peak 2:	0.000	0.0	Z-Average (d.nm): 653.4 PDI: 0.086
Peak 3:	0.000	0.0	Count Rate (kcps): 290.7

Note: Z-average and PDI are based on the intensity distribution



General Notes: Average result created from record number(s): 1 2 3 4 5 6

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	8
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator:
Malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021 BRJ

Approved:

28 JULI 2021

WY

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32030
Record Number: 9
27 jul 2021 14:12



CERTIFICATE OF ANALYSIS

Customer: CphNano
Material tested: Polystyrene Latex beads 1100nm
Batch: 2021/07/23
Internal number: 32031

Analytical technique: Dynamic Light Scattering
Method of analysis: Analyseplan
Internal quality level: GMP

	D _{10%} (nm)**	D _{50%} (nm)**	D _{90%} (nm)**	Z-average*	PDI*
Run 1	527	873	1250	789.6	0.147
Run 2	736	953	1230	873.5	0.04
Average	620	922	1240	831.5	0.093
Specification	-	-	-	-	-
Evaluation	-	-	-	-	-

*The parameter is based on the intensity size distribution

**The parameter is based on the volume size distribution

Written by: Bent Duusgaard

Date: 28 JULI 2021

Reviewed by: Wenbo Wang

Date: 28 JULI 2021

Approved by: Wenbo Wang
QC

Date: 28 JULI 2021

The validity of the method is the responsibility of the sponsor
Quality agreement not in place

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32031 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:32:43

Userid: brj

	D10%(V): 849 nm	D50%(V): 1100 nm	D90%(V): 1420 nm			
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and Pdl are based on the intensity distribution
Peak 1:	1112	100.0	194.5			
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	272.1	



General Notes: 32031 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorption:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	9
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

28 JULI 2021

Operator: 27 JULI 2021 BRJ
malvern Panalytical
www.malvernpanalytical.com

Approved: BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32031
Record Number: 27 jul 2021 14:12

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32031 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:34:57

Userid: brj

	D10%(V): 430 nm	D50%(V): 609 nm	D90%(V): 841 nm
	Diam. (nm)	% Volume	Width (nm)
Peak 1:	623.1	100.0	149.4
Peak 2:	0.000	0.0	0.000
Peak 3:	0.000	0.0	0.000
		Z-Average (d.nm):	592.9
		PdI:	0.264
		Count Rate (kcps):	269.5

Size Distribution by Volume



Record 2: 32031 Batch 2021/07/23

General Notes: 32031 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	9
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

28 JULI 2021

Operator: brj

27 JULI 2021 brj

Approved: brj

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32031
Record Number: 2
27 jul 2021 14:12



Sample Name: 32031 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27.juli 2021 11:37:10

Userid: brj

	D10%(V): 677 nm	D50%(V): 897 nm	D90%(V): 1170 nm
Peak 1:	Diam. (nm) 905.8	% Volume 100.0	Width (nm) 167.7
Peak 2:	0.000	0.0	Z-Average (d.nm): 800.2 PDI: 0.071
Peak 3:	0.000	0.0	Count Rate (kcps): 271.1

Note: Z-average and PDI are based on the intensity distribution



General Notes: 32031 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	9
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021 BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved:

28 JULI 2021 WJP

File name: 32031
Record Number: 3
27 jul 2021 14:12

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32031 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:43:54

Userid: brj

D10%(V): 721 nm			D50%(V): 924 nm			D90%(V): 1190 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	859.6	Note: Z-average and PDI are based on the intensity distribution	
Peak 1:	933.9	100.0	158.5	PDI:	0.088		
Peak 2:	0.000	0.0	0.000				
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		263.8	



General Notes: 32031 Run 2

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000 nm	
			Number of measurements:	3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	9
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	13

Operator:

27 JULI 2021 *BRJ*

Approved:

28 JULI 2021 *WJP*

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32031
Record Number: 4
27 jul 2021 14:12

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32031 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:46:07

Userid: brj

	D10%(V): 762 nm	D50%(V): 999 nm	D90%(V): 1280 nm			
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution
Peak 1:	1010	100.0	182.1	904.9	0.012	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 265.3		



General Notes: 32031 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	9
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator:
alvern Panalytical
www.alvernpanalytical.com

27 JULI 2021 *BRJ*

Approved:

28 JULI 2021 *WJ*

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32031
Record Number: 5
27 jul 2021 14:12

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32031 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:48:20

Userid: brj

	D10%(V): 732 nm	D50%(V): 944 nm	D90%(V): 1210 nm
Peak 1:	Diam. (nm) 954.1	% Volume 100.0	Width (nm) 162.4
Peak 2:	0.000	0.0	Z-Average (d.nm): 856.0 PDI: 0.019
Peak 3:	0.000	0.0	Count Rate (kcps): 264.1

Note: Z-average and PDI are based on the intensity distribution



General Notes: 32031 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	9
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved:

28 JULI 2021

WJ

File name: 32031
Record Number: f
27 jul 2021 14:12

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32031 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:48:41

Userid: brj

D10%(V): 527 nm			D50%(V): 873 nm	D90%(V): 1250 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	789.6	Note: Z-average and Pdl are based on the intensity distribution
Peak 1:	880.3	100.0	263.8	PdI:	0.147	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):	272.1	



General Notes: Average result created from record number(s): 1 2 3

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	9
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021 BRJ

Approved:

28 JULI 2021

W.M.

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32031
Record Number: 7
27 jul 2021 14:12

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32031 Batch 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:48:50

Userid: brj

	D10%(V): 736 nm	D50%(V): 953 nm	D90%(V): 1230 nm			
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution
Peak 1:	965.9	100.0	171.0	873.5	0.040	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000			
				Count Rate (kcps):	263.8	



General Notes: Average result created from record number(s): 4 5 6

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	9
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	13

Operator: 27 JULI 2021 brj
alvern Panalytical
www.alvernpanalytical.com

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved: 28 JULI 2021 W.M.

File name: 32031
Record Number: 8
27 jul 2021 14:12

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS

Sample Name: 32031 Batch 2021/07/23


SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 11:49:02

Userid: brj

	D10%(V): 620 nm	D50%(V): 922 nm	D90%(V): 1240 nm
Peak 1:	Diam. (nm) 923.1	% Volume 100.0	Width (nm) 226.4
Peak 2:	0.000	0.0	Z-Average (d.nm): 831.5 PDI: 0.093
Peak 3:	0.000	0.0	Count Rate (kcps): 272.1

Note: Z-average and PDI are based on the intensity distribution



General Notes: Average result created from record number(s): 1 2 3 4 5 6

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
Auto Position Enabled:	False	Size Measure Delay (s):	0
Auto Attenuate Enabled	True	Measurement Position (mm):	4.65
Auto Size Measurement Time:	True	Attenuator:	9
		Duration (s):	10 Extend duration for large particles: False
		Size Runs:	13

Operator: 27 JULI 2021 BRJ

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved: 28 JULI 2021

BRJ

File name: 32031
Record Number: 9
27 jul 2021 14:12



CERTIFICATE OF ANALYSIS

Customer: CphNano
Material tested: Polystyrene Latex beads 3000nm
Batch: 2021/07/23
Internal number: 32032

Analytical technique: Dynamic Light Scattering
Method of analysis: Analyseplan
Internal quality level: GMP

	D _{10%} (nm)**	D _{50%} (nm)**	D _{90%} (nm)**	Z-average*	PDI*
Run 1	1040	1690	4030	2698	0.365
Run 2	1150	1480	2290	2756	0.334
Average	1110	1550	3570	2727	0.349
Specification	-	-	-	-	-
Evaluation	-	-	-	-	-

*The parameter is based on the intensity size distribution

**The parameter is based on the volume size distribution

Written by: Wenbo Wang

Date: 28 JULI 2021

Reviewed by: Wenbo Wang

Date: 28 JULI 2021

Approved by: Wenbo Wang
QC

Date: 28 JULI 2021

The validity of the method is the responsibility of the sponsor
Quality agreement not in place

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27.juli 2021 12:47:28

Userid: brj

D10%(V): 997 nm			D50%(V): 1260 nm			D90%(V): 1600 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution	
Peak 1:	1272	100.0	196.0			2751	
Peak 2:	0.000	0.0	0.000				
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		220.2	



General Notes: 32032 Run 1

Cell Description:	Disposable sizing cuvette			
Material RI:	1.59		Dispersant Name:	Water
Material Absorbtion:	0.010		Dispersant RI:	1.330
Analysis Model:	General Purpose		Viscosity (cP):	0.8872
Lower Size Threshold:	0.050		Temperature (°C):	25.0
Upper Size Threshold:	0.010		Equilibration Time Set (min):	180
Size range:	0.6000	to	6000 nm	
			Number of measurements:	3
			Size Measure Delay (s):	0
Auto Position Enabled:	False		Measurement Position (mm):	4.65
Auto Attenuate Enabled	True		Attenuator:	6
Auto Size Measurement Time:	True		Duration (s):	10
			Extend duration for large particles:	False
			Size Runs:	14

Operator: 27 JULI 2021 BRJ

Approved:

28 JULI 2021

WJ

File name: 32032
Record Number:
27 jul 2021 14:13

Zetasizer Ver. 7.13
Serial Number : MAL500686

Size Distribution Report by Volume




Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 12:49:52

Userid: brj

D10%(V): 1410 nm			D50%(V): 2100 nm			D90%(V): 3490 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution	
Peak 1:	2289	100.0	792.3			2779	
Peak 2:	0.000	0.0	0.000				
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		211.2	



General Notes: 32032 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	14

Operator:
 alvern Panalytical
 www.alvernpanalytical.com

27 JULI 2021 BRJ

Approved:

 Zetasizer Ver. 7.13
 Serial Number : MAL500686

28 JULI 2021 WJ

 File name: 32032
 Record Number: 2
 27 jul 2021 14:13

Size Distribution Report by Volume



**PARTICLE
ANALYTICAL**
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 12:52:15

Userid: brj

	D10%(V): 866 nm	D50%(V): 2760 nm	D90%(V): 5200 nm
	Diam. (nm)	% Volume	Width (nm)
Peak 1:	825.6	15.3	145.0
Peak 2:	1655	20.2	275.6
Peak 3:	3683	64.5	1207
		Z-Average (d.nm):	2563
		PDI:	0.238
		Count Rate (kcps):	191.0



Record 3: 32032 Baych 2021/07/23

General Notes: 32032 Run 1

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	14

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021

Approved:

28 JULI 2021

WJP

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32032
Record Number: 3
27 jul 2021 14:13

Size Distribution Report by Volume




Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 12:59:22

Userid: brj

	D10%(V): 1160 nm	D50%(V): 1430 nm	D90%(V): 1750 nm
	Diam. (nm)	% Volume	Width (nm)
Peak 1:	1443	100.0	192.2
Peak 2:	0.000	0.0	0.000
Peak 3:	0.000	0.0	0.000
	Z-Average (d.nm):	3136	Note: Z-average and PDI are based on the intensity distribution
	PDI:	0.444	
	Count Rate (kcps):	217.0	



General Notes: 32032 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	15

Operator:

27 JULI 2021

Approved:

28 JULI 2021

verj

File name: 32032
Record Number: 4
27 jul 2021 14:13

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 13:01:56

Userid: brj

D10%(V): 1230 nm			D50%(V): 1570 nm			D90%(V): 2270 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and Pdl are based on the intensity distribution		
Peak 1:	1668	100.0	427.5			2690		
Peak 2:	0.000	0.0	0.000			0.265		
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		191.6		



General Notes: 32032 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	15

Operator:
alvern Panalytical
www.malvernpanalytical.com

27 JULI 2021 BRJ

Approved:

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32032
Record Number: f
27 jul 2021 14:13

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 13:04:29

Userid: brj

D10%(V): 1060 nm			D50%(V): 1460 nm			D90%(V): 4910 nm		
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution		
Peak 1:	1425	81.6	357.4					
Peak 2:	4966	18.4	899.5					
Peak 3:	0.000	0.0	0.000	Count Rate (kcps): 173.8				



General Notes: 32032 Run 2

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorption:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	15

Operator:
alvern Panalytical
www.alvernpanalytical.com

27 JULI 2021 *BRJ*

Zetasizer Ver. 7.13
Serial Number : MAL500686

Approved:

28 JULI 2021 *W.M.*

File name: 32032
Record Number: f
27 jul 2021 14:13




Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 13:06:01

Userid: brj

D10%(V): 1040 nm			D50%(V): 1690 nm			D90%(V): 4030 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and PDI are based on the intensity distribution	
Peak 1:	1388	59.0	332.3				
Peak 2:	3243	41.0	1098				
Peak 3:	0.000	0.0	0.000	Count Rate (kcps):		220.2	



General Notes: Average result created from record number(s): 1 2 3

Cell Description: Disposable sizing cuvette

Material RI: 1.59

Dispersant Name: Water

Material Absorbtion: 0.010

Dispersant RI: 1.330

Analysis Model: General Purpose

Viscosity (cP): 0.8872

Lower Size Threshold: 0.050

Temperature (°C): 25.0

Upper Size Threshold: 0.010

Equilibration Time Set (min): 180

Size range: 0.6000 to 6000 nm

Number of measurements: 3

Size Measure Delay (s): 0

Auto Position Enabled: False

Measurement Position (mm): 4.65

Auto Attenuate Enabled: True

Attenuator: 6

Auto Size Measurement Time: True

Duration (s): 10 Extend duration for large particles: False

Size Runs: 14

Operator: 27 JULI 2021 BRJ

28 JULI 2021

Approved:

 File name: 32032
 Record Number: 7
 27 jul 2021 14:13

 Zetasizer Ver. 7.13
 Serial Number : MAL500686

Size Distribution Report by Volume



PARTICLE
ANALYTICAL
EXPERTS IN SIZE AND CRYSTALS


Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 13:06:11

Userid: brj

D10%(V): 1150 nm			D50%(V): 1480 nm		D90%(V): 2290 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PDI:	Note: Z-average and Pdl are based on the intensity distribution
Peak 1:	1527	94.0	383.5			
Peak 2:	5020	6.0	836.9			
Peak 3:	0.000	0.0	0.000			
				Count Rate (kcps):	217.0	



General Notes: Average result created from record number(s): 4 5 6

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10 Extend duration for large particles: False
		Size Runs:	15

Operator:
malvern Panalytical
www.malvernpanalytical.com

27 JULI 2021 BRJ

Approved:

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686

File name: 32032
Record Number: 8
27 jul 2021 14:13

Size Distribution Report by Volume




Sample Name: 32032 Baych 2021/07/23

SOP Name: OQ - 173 degree, 633nm laser, 60 nm dilute latex.sop

Measurement Date and Time: 27. juli 2021 13:06:47

Userid: brj

D10%(V): 1110 nm			D50%(V): 1550 nm		D90%(V): 3570 nm	
	Diam. (nm)	% Volume	Width (nm)	Z-Average (d.nm):	PdI:	Note: Z-average and PdI are based on the intensity distribution
Peak 1:	1943	100.0	1103	2727	0.349	
Peak 2:	0.000	0.0	0.000			
Peak 3:	0.000	0.0	0.000		Count Rate (kcps): 220.2	



General Notes: Average result created from record number(s): 1 2 3 4 5 6

Cell Description:	Disposable sizing cuvette		
Material RI:	1.59	Dispersant Name:	Water
Material Absorbtion:	0.010	Dispersant RI:	1.330
Analysis Model:	General Purpose	Viscosity (cP):	0.8872
Lower Size Threshold:	0.050	Temperature (°C):	25.0
Upper Size Threshold:	0.010	Equilibration Time Set (min):	180
Size range:	0.6000 to 6000 nm	Number of measurements:	3
		Size Measure Delay (s):	0
Auto Position Enabled:	False	Measurement Position (mm):	4.65
Auto Attenuate Enabled	True	Attenuator:	6
Auto Size Measurement Time:	True	Duration (s):	10
		Extend duration for large particles:	False
		Size Runs:	14

Operator: 27 JULI 2021 BRJ

Approved:

File name: 32032
Record Number: 9
27 jul 2021 14:13

28 JULI 2021

Zetasizer Ver. 7.13
Serial Number : MAL500686