

Access the Breadth of Sample Preparation Protocols to Awaken True Functional Immune Biology

Precedent Protocols and Processes ensure that enrichment, stimulation and labeling will result in maximizing the correlative data achieved, while minimizing non-specific impact on the cells themselves.


Highlights:

- Published enrichment, stimulation, and staining protocols throughout immunotherapy types and indications help achieve your goals, in multiple publications and presentations
- Protocols that minimize impacts on the live cells that are analyzed on Bruker's IsoLight and IsoSpark platform
- Guidelines for most commonly used cell types and applications



Protocols to engage polyfunctional T cells in checkpoint and combination therapies

PSI correlative biomarkers in a wide body of work



Proven protocols awaken correlative biology



From accessible blood sample, requiring lower volume of cells

Figure 1 | Our polyfunctional biomarkers reveal the uniquely correlative T-cell specific mechanism in TILS/Bone Marrow microenvironment in a variety of therapeutic areas: checkpoint therapy, combination agonists, Double/ triple combination checkpoint, and Immune fitness. PSI reveals these unique correlates in the blood as well, which presents a more actionable & accessible source of sample in: Cellular cancer vaccines, Checkpoint therapy, Combination vaccines, Bi-specifics, and TCR-T cell Therapies. Our proven Protocols Awaken Correlative Cytokine Biology with stimulation to reveal true functional t cell biology..

Protocol frameworks where Bruker systems achieved published correlates: workflows

Checkpoint & Combination Therapies						
IsoCode Single Cell Protocols	Immune Cell Type	Immunotherapy Type	Indication	Enrichment	Stimulation Type & Time	Polyfunctional Correlate/ Difference
PBMC: T Cell Protocol (H)	Human PBMC ¹³	Pre and post cancer vaccine ¹³	Cancer: Pancreatic Cancer	Miltenyi Beads	CD3 / CD28: 24 hours	Survival correlate
PBMC: T Cell Protocol (H)	Human PBMC ¹¹	Pre and post TLR Agonist Therapy ¹¹	Cancer: Sarcoma	Miltenyi Beads	CD3/CD28: 7 hours	TCR Seq Correlate
PBMC: T Cell Protocol (H)	Humanized Mouse PBMC ²¹	Nanoparticle cancer vaccine ¹³	Cancer: Solid Tumor	Miltenyi Beads	CD3 / CD28: 24 hours	Mouse In Vivo Correlate
T Cell Protocol (H)	Human Bone Marrow ¹⁹	Combination Checkpoint Therapy	Cancer: Leukemia	Miltenyi Beads	CD3 / CD28: 24 hours	Response Correlates
T Cell Protocol (H)	Human Bone Marrow ¹⁹	Combination Checkpoint Therapy	Cancer: Leukemia	Miltenyi Beads	CD3 / CD28: 24 hours	Newly Diagnosed & Relapse
TILs Protocol (H)	Human TILs ¹²	Post Checkpoint Therapy ¹²	Cancer: Solid Tumor	Miltenyi Beads	CD3 / CD28: 24 hours	Responder correlate
TILs Protocol (H)	Human TILs ²⁰	Triple Checkpoint Therapy	Cancer: Ovarian Cancer	Miltenyi Beads	CD3 / CD28: 24 hours	Functionality Correlates
Bispecific Protocol (H)	Human PBMC ¹⁸	Bispecific T cell engager	Cancer: Solid Tumor	Miltenyi Beads	Bispecific T cell engagers: 36 hours	Preclinical Correlate
TILs Protocol (M)	Mouse TILs ¹⁶	Post combination therapy	Cancer: Melanoma	Miltenyi Beads	CD3/CD28: 24 hours	Mouse In vivo response

Figure 2 | The goal for starting with a straightforward protocol that is recommended by Bruker is to ensure success and drive optimization. These protocols have been highly successful in detecting polyfunctional differences in the past (see Technology Note: Validation of Sample Preparation). We recommend not deviating from published protocols, but at the same time simplifying certain aspects that should not impact the success of the study, to ensure success for first time IsoLight users.

Protocols to engage polyfunctional T cells and immune cells in cell therapies

PSI correlative biomarkers in a wide body of work



- CAR T Therapies
- TCR-T Therapies
- NK cell Therapies
- Bispecific Therapies
- CRISPR edited Therapies



- Checkpoint therapy
- Combination vaccines
- Bi-specific CAR T optimization
- TCR-T cell Therapies

Proven protocols awaken correlative biology



Directly from pre-infusion product, or blood, requiring lower volume

Figure 3 | Uncover True functional t cell biology and discover the subsets of cells that make the difference in engineered cell therapies: CAR T , TCR T, NK Cell, Bispecific, and CRISPR edited therapies. PSI reveals these unique correlates in Immun-oncology: Checkpoint therapy, combination vaccines, BiSpecific CAR T Optimization and TCR T cell therapies. Our proven Protocols Awaken Correlative Cytokine Biology with stimulation to reveal true functional t cell biology directly from the preinfusion product or patient blood samples.

Protocol frameworks where Bruker systems achieved published correlates: workflows

Cell Engineering & Therapy						
IsoCode Single Cell Protocols	Immune Cell Type	Immunotherapy Type	Indication	Enrichment	Stimulation Type & Time	Polyfunctional Correlate/ Difference
CAR-T Protocol (H)	Human CD19 CAR- T Cells ²	CAR-T Therapy	Cancer: Non- hodgkins Lymphoma	Miltenyi Beads	CD19 K562 target cell: 20 hours	Responder Correlate
CAR-T Protocol (H)	Human CD19 CAR- T Cells ⁶	CAR-T Therapy	Cancer: B cell ALL	Miltenyi Beads	CD19 Beads: 24 hours	Product Heterogeneity Differences
TCR-T Protocol (H)	Human TCR- Engineered Mart-1 T cells ⁷	TCR-T Therapy	Cancer: Melanoma	Tetramer	MART-1 tetramer: 12 hours	Tumor relapse correlate
NK Cell Protocol (H)	Human NK Cells ¹⁰	NK Cell Therapy	Cancer: Leukemia	Miltenyi Beads	IL-15 O/N + IL-12 / IL-18	Mouse tumor regression correlate
Bispecific Protocol (H)	Human Bispecific CAR-T Cells ¹⁸	Bispecific CAR-T Therapy	Cancer: B cell ALL	Miltenyi Beads	CD19/CD22 K562: 20 hours	Process Optimization Characterization
TILs Protocol (M)	Mouse TILs ¹⁶	Combination ACT	Cancer: Melanoma	Miltenyi Beads	CD3/CD28: 24 hours	Mouse In vivo response
CAR-T Protocol (M)	Mouse CD19 CAR- T Cells ¹⁷	CAR-T Therapy	Cancer: B-cell acute lymphoblastic leukemia	Miltenyi Beads	3T3-mCD19 Target cell: 24 hours	Young vs Aging Functional Correlate

Figure 4 | The goal for starting with a straightforward protocol that is recommended by Bruker is to ensure success and drive optimization. These protocols have been highly successful in detecting polyfunctional differences in the past (see Technology Note: Validation of Sample Preparation). We recommend not deviating from published protocols, but at the same time simplifying certain aspects that should not impact the success of the study, to ensure success for first time IsoLight users.

Protocol frameworks where Bruker systems achieved published correlates: workflows

Vaccine						
IsoCode Single Cell Protocols	Immune Cell Type	Therapy Type	Indication	Enrichment	Stimulation Type & Time	Polyfunctional Correlate/ Difference
PBMC: T Cell Protocol (H)	Human PBMC / T cells ¹³	Cancer Vaccine	Cancer: Pancreatic Cancer	Miltenyi Beads	CD3/CD28: 24 hours	Survival correlates
PBMC: T Cell Protocol (H)	Humanized Mouse T Cells ²¹	Cancer Vaccine	Cancer: Melanoma	Miltenyi Beads	CD3/CD28: 24 hours	In Vivo Response
PBMC: T Cell Protocol (M)	Mouse T cell ²²	Malaria Vaccine	Infectious Disease: Malaria	Miltenyi Beads	CD3/CD28: 24 hours	In Vivo Response
Inflammation & Autoimmunity						
IsoCode Single Cell Protocols	Immune Cell Type	Therapy Type	Indication	Enrichment	Stimulation Type & Time	Polyfunctional Correlate/ Difference
Monocyte Protocol (H)	Human Monocyte in healthy and diseased patients ¹⁴		Autoimmune & CNS: Multiple Sclerosis	Miltenyi Beads	PC3 or LPS: 24 hours	MS correlate
CAR-T Protocol (H)	Human CAR T cell Product ²	CAR-T Therapy	Immune Related Adverse Events: Neurotoxicity	Miltenyi Beads	CD19-K562: 20 hours	Grade 3+ CRS/ Neurotoxicity
T Cell Protocol (H)	Human T cells ²³		Autoimmune: Systemic Lupus Erythematosus	Flow Cytometry	No Stimulation	Inflammatory state correlate
NK Cell Protocol (H)	Human NK cells, no therapy ⁸		Autoimmune: IBD and Crohn's disease	Flow Cytometry	PMA/Ionomycin: 12 hours	Inflammatory state correlate

Figure 5 | The goal for starting with a straightforward protocol that is recommended by Bruker is to ensure success and drive optimization. These protocols have been highly successful in detecting polyfunctional differences in the past (see Technology Note: Validation of Sample Preparation). We recommend not deviating from published protocols, but at the same time simplifying certain aspects that should not impact the success of the study, to ensure success for first time IsoLight users.

Bruker Protocol Frameworks to Provide Users with Precedent Workflows

Bruker systems have been used in a variety of research areas and contexts, as seen in Figure 2, 4, 5 to achieve understanding of correlates to in vivo activity in patients and mice. The below figure provides an overview of the variety of research areas, cell types, and indications where Bruker systems have been used to publish data to achieve correlates to response, i.e. patient differences, that explain mechanism in vivo. The Stimulation type and enrichment

type can act as a guide, with cited publications and presentations, that users can reference.

See Figure 2, 4, 5 for a variety of contexts in which our protocols for stimulation, enrichment, and labeling were used to achieve correlative outcome.

Conclusion

- Published enrichment, stimulation, and staining protocols throughout immunotherapy types and indications can help provide a framework to achieve goals of obtaining single-cell data, in multiple publications and presentations
- These enrichment and stimulation protocols minimize impacts on the live cells that are analyzed in Bruker's systems
- The type of sample preparation protocols and strategies have helped researchers achieve critical correlates to in vivo data in the past as well
- While every user project is specific in terms of requirements, (e.g., stimulation time), precedent datasets can act as a rubric for choice of protocol.

References

- Lu Y, Xue Q, Eisele MR, Sulistijo ES, Brower K, Han L, Amir el AD, Pe'er D, Miller-Jensen K, Fan R (2015) Highly multiplexed profiling of single-cell effector functions reveals deep functional heterogeneity in response to pathogenic ligands. *Proc Natl Acad Sci U S A* 112 (7):E607-615. doi:10.1073/pnas.1416756112
- Rossi J, Paczkowski P, Shen YW, Morse K, Flynn B, Kaiser A, Ng C, Gallatin K, Cain T, Fan R, Mackay S, Heath JR, Rosenberg SA, Kochenderfer JN, Zhou J, Bot A (2018) Preinfusion polyfunctional anti-CD19 chimeric antigen receptor T cells are associated with clinical outcomes in NHL. *Blood* 132 (8):804-814. doi:10.1182/blood-2018-01-828343
- Precopio ML, Betts MR, Parrino J, Price DA, Gostick E, Ambrozak DR, Asher TE, Douek DC, Harari A, Pantaleo G, Bailer R, Graham BS, Roederer M, Koup RA (2007) Immunization with vaccinia virus induces polyfunctional and phenotypically distinctive CD8(+) T cell responses. *J Exp Med* 204 (6):1405-1416. doi:10.1084/jem.20062363
- Ding ZC, Huang L, Blazar BR, Yagita H, Mellor AL, Munn DH, Zhou G (2012) Polyfunctional CD4(+) T cells are essential for eradicating advanced B-cell lymphoma after chemotherapy. *Blood* 120 (11):2229-2239. doi:10.1182/blood-2011-12-398321
- Darrah PA, Patel DT, De Luca PM, Lindsay RW, Davey DF, Flynn BJ, Hoff ST, Andersen P, Reed SG, Morris SL, Roederer M, Seder RA (2007) Multifunctional TH1 cells define a correlate of vaccine-mediated protection against Leishmania major. *Nat Med* 13 (7):843-850. doi:10.1038/nm1592
- Xue Q, Bettini E, Paczkowski P, Ng C, Kaiser A, McConnell T, Kodrasi O, Quigley MF, Heath J, Fan R, Mackay S, Dudley ME, Kassim SH, Zhou J (2017) Single-cell multiplexed cytokine profiling of CD19 CAR-T cells reveals a diverse landscape of polyfunctional antigen-specific response. *J Immunother Cancer* 5 (1):85. doi:10.1186/s40425-017-0293-7
- Ma C, Cheung AF, Chodon T, Koya RC, Wu Z, Ng C, Avramis E, Cochran AJ, Witte ON, Baltimore D, Chmielewski B, Economou JS, Comin-Anduix B, Ribas A, Heath JR (2013) Multifunctional T-cell analyses to study response and progression in adoptive cell transfer immunotherapy. *Cancer Discov* 3 (4):418-429. doi:10.1158/2159-8290.CD-12-0383
- Lin L, Ma C, Wei B, Aziz N, Rajalingam R, Yusung S, Erlich HA, Trachtenberg EA, Targan SR, McGovern DP, Heath JR, Braun J (2014) Human NK cells licensed by killer Ig receptor genes have an altered cytokine program that modifies CD4+ T cell function. *J Immunol* 193 (2):940-949. doi:10.4049/jimmunol.1400093
- Fujiwara M, Anstadt EJ, Flynn B, Morse K, Ng C, Paczkowski P, Zhou J, Mackay S, Wasko N, Nichols F, Clark RB (2018) Enhanced TLR2 responses in multiple sclerosis. *Clin Exp Immunol* 193 (3):313-326. doi:10.1111/cei.13150
- Zhu H, Blum RH, Wu Z, Bahena A, Hoel HJ, Ask EH, Guan KL, Malmberg KJ, Kaufman DS. Notch activation rescues exhaustion in CISH-deleted natural killer cells to promote in vivo persistence and enhance anti-tumor activity. *ASH Annual Meeting* 2018.
- Seo YD, Zhou J, Morse K, Patino J, Mackay S, Kim EY, Conrad III EU, O'Malley RB, Cranmer L, Lu H, Hsu FJ, Xu Y, Loggers E, Hain T, Pillarisetty VG, Kane G, Riddell S, Meulen J, Jones RL, Pollack SM. Intratumoral (IT) injection of the toll-like receptor 4 (TLR4) agonist G100 induces a clinical response and a T cell response locally and systemically. *Journal of Clinical Oncology*, 36, Suppl 5S (2018).
- Mackay S, Flynn B, Morse K, Paczkowski P, Chen J, Liu D, Bacchiocchi A, Heath JR, Fan R, Sznol M, Halaban R, Zhou J. Single-cell PSI of CD8+ TILs in melanoma shows uniquely sensitive correlates with response to anti-PD-1/CTLA4 therapy, where histology and serum cytokines were unable to detect significant associations. *SITC Annual Meeting* 2018.
- Mackay S, Flynn B, Chen J, Paczkowski P, Jaffee E, Zheng L, and Zhou J. Single-Cell Polyfunctionality of CD4+ T Cells Shows Promise as a Predictor of Overall Survival of Pancreatic Cancer Patients Treated with GVAX Vaccine. *FOCIS* 2018.
- Fujiwara M, Anstadt EJ, Flynn B, Morse K, Ng C, Paczkowski P, Zhou J, Mackay S, Wasko N, Nichols F, Clark RB. Enhanced TLR2 Responses in Multiple Sclerosis. *Clinical and Experimental Immunology* (2018).
- Lin L, Ma C, Wei B, Aziz N, Rajalingam R, Yusung S, Erlich HA, Trachtenberg EA, Targan SR, McGovern DPB, Heath JR, and Braun J. Human NK Cells Licensed by Killer Ig Receptor Genes have an Altered Cytokine Program that Modifies CD4+ T Cell Function. *Journal of Immunology*, 193, 940-9 (2014).
- Parisi G, Saco J, Salazar F, Krystofinski P, Tsoi J, Zhang R, Puig Saus C, Zhou J, Hu-Lieskova S, Comin-Anduix B, Wu A, Charych DH, Ribas A. Enhanced Expansion and Tumor Targeting of Adoptively Transferred T Cells with NKTR-214. Session PO.CL06.03 - Adoptive Cell Therapy 3. American Association for Cancer Research Annual Meeting 2018.
- Kotani H, Li Gongbo, Yao Jiqiang, Mesa TE, Chen J, Boucher JC, Yoder SJ, Zhou J, Davila ML. Aged CAR T cells exhibit enhanced cytotoxicity and effector function but shorter persistence and less memory-like phenotypes. *ASH Annual Meeting* 2018
- Mackay S, Paczkowski P, Flynn B, Morse K, Liu D, Coupet TA, Godbersen C, Sentman CL, Zhou J. Single-cell proteomic analysis of T cells stimulated by Bi-specific T cell Engagers shows a robust and unique polyfunctional secretion profile. *SITC Annual Meeting* 2018..
- Daver N, Kantarjian H, Mackay S, Flynn B, Basu S, Garcia-Manero G, Alfayez M, Konopleva M, Matthews J, Kornblau S, Jabbour E, Zhou J, Andreeff M. Polyfunctionality determined by single-cell proteomics of bone marrow-derived CD4 T cells from patients with acute myeloid leukemia identifies patients responding to anti-PD-1-based therapy and discovers profound T cell defect in mutant TP53 disease. *AACR Conference* 2019.
- Kaufmann JK, Flynn B, Morse K, Speranza MC, Zhou J, Ramaswamy S, Mackay S, Coleman KG. Triple checkpoint blockade targeting PD-1, TIM-3, and LAG-3 reinvigorates ovarian cancer-infiltrating T cells by increasing T cell polyfunctionality and effector function. *AACR Conference* 2019.
- Mackay S, Morse K, Paczkowski P, Huang J, Tsuji M, and Zhou J. Single-Cell Highly Multiplexed Proteomics Identifies Novel Polyfunctional Human CD8+ T Cell Signatures Induced by a Nanoparticle-Based Melanoma Vaccine in Human Immune System Mice. *FOCIS* 2018.
- Zhou J, Kaiser A, Ng C, Karcher R, McConnell T, Paczkowski P, Fernandez C, Zhang M, Mackay S, and Tsuji M. CD8+ T-Cell Mediated Anti-Malaria Protection Induced by Malaria Vaccines; Assessment of Hepatic CD8+ T Cells by SCBC Assay. *Human Vaccines & Immunotherapeutics*, 13 (7), 1625-1629 (2017)
- Fan R et al. Phenotypic and Functional Heterogeneity of cTfh Cells in Systemic Lupus Erythematosus. *FOCIS* 2016.
- Fan R. Single Cell Analysis of Cytokine Production. *ASH* 2016.