

# COVID-19 Highlight: Running Plasma and CSF Samples on the CodePlex Secretome Solution

The following workflow has been used in generating highly multiplexed population cytokine data and serves as a guideline for running COVID-19 patient plasma and cerebrospinal fluid (CSF) samples using the CodePlex Secretome solution.

---

## **In this Technical Note we outline:**

- Safety precautions
- Preparation of patient plasma and CSF samples for loading



### Proper Safety Precautions


- A disposable laboratory coat should always be worn when handling COVID19 samples.
  - A secondary coat can be worn for extra precaution. Coat can be worn in the reverse manner as the first.
- An N-95 mask should always be worn when handling COVID19 samples.
- A face shield or safety goggles should be worn when handling COVID19 samples.
- Disposable sleeves should be worn when handling COVID19 samples.
  - These sleeves should go under the disposable gown/s but above disposable gloves.
  - Sleeves and gloves should be taped together to prevent skin exposure.
- Disposable nitrile gloves should always be worn when handling COVID19 samples.
  - Gloves can be double gloved barring dexterity is not compromised.
- Remove all disposable PPE and place into larger biohazard container, seal biohazard container.
- All additional lab safety procedures should be followed according to each lab's requirements dependent on sample type.

## Prep, Run, Analyze

### Preparation of plasma and CSF samples for loading

**Note:** All steps should be done in a biological safety hood.

#### Codeplex workflow for patient plasma and CSF



Collected Plasma and CSF samples are thawed and loaded on the CodePlex Secretome chip. The IsoLight system automates the entire ELISA workflow.

- 1. Thaw vials of CSF.** Allow CSF to thaw and come to a complete liquid form (30-60 minutes) at room temperature.
- 2. Thaw vials of plasma.** Allow plasma to thaw and come to a complete liquid form (30-60 minutes) at room temperature.
- 3. Thaw CodePlex chips.** Allow CodePlex chips to thaw for 30-60 minutes prior to media loading at room temperature.
- 4. Pipet up and down gently.** Place excess sample vials in -80°C freezer.
- 5. Load CSF onto CodePlex chip.** CSF IS NOT TO BE DILUTED. Background wells are to be filled with DI water.
- 6. Load plasma onto second CodePlex chip.**
- 7. Perform a 1:1 dilution of plasma and 2% BSA/PBS.** Mix gently and load on to third CodePlex chip.
- 8. Perform a 1:2 dilution of plasma and 2% BSA/PBS.** Mix gently and load on to fourth CodePlex chip.
- 9. Load all chips onto IsoLight and run experiment.** Load CodePlex reagents prior to insertion of CodePlex chips into the IsoLight. Before starting the run, add 100mL of common undiluted bleach into the IsoLight waste container.
- 10. Dispose of all hazardous materials.** After the run, remove CodePlex chips and place into a sharps biohazard container, seal container and place sharps container into a larger biohazard container.
- 11. Run an 8 chip cleaning script on the IsoLight.** Run a secondary 8 chip cleaning script as precaution. Upon completion of two cleaning runs, remove waste container from IsoLight and top off with common undiluted bleach in a biological safety hood and allow to sit for 2-3 hours.

**ALL ADDITIONAL LAB SAFETY PROCEDURES SHOULD BE FOLLOWED ACCORDING TO EACH LAB'S REQUIREMENTS DEPENDENT ON SAMPLE TYPE.**

**DISPOSE OF WASTE CONTAINER INTO A BIOHAZARD CONTAINER.**