

Proning Therapy Ventilated and Non-Ventilated Orders

MEDITECH 5.6

Proning Therapy Ventilated

Note: The Proning Therapy Order can be set up to be ordered as a series, with the expectation of modifying the order if a different proning duration is desired.

Discontinue the order if therapy is not desired for that day. This should not function as a PRN order.

Screen Sample

Procedure Ordered
Proning Therapy Ventilated

Pri	Service Date/Time	Series	Directions	Qty
	07/20			

Proning duration:

1	8 hours	5	20 hours	
2	12 hours	6	Other	
3	16 hours			
4	18 hours			

Proning duration: *

Proning other duration:

Proning comment:

Ok Cancel Help Prev Next

Proning Therapy Non-Ventilated

Screen Sample

Procedure Ordered
Proning Therapy Non-Ventilated

Pri	Service Date/Time	Series	Directions	Qty
	07/20			

Proning duration:

1	1 hour	5	5 hour	9 Other
2	2 hour	6	6 hour	
3	3 hour	7	7 hour	
4	4 hour	8	8 hour	

Proning duration: *

Proning other duration:

Proning comment:

Ok Cancel Help Prev Next