

Standardized Central Line Dressing Kit (a collaboration between CSG

A new kit designed by nursing, for nurses, is now available.

- ⇒ Two step kit contained in one package
- ⇒ Step 1 for removal of the current dressing
- ⇒ Step 2 for application of the new dressing
- ⇒ Eliminates waste with all necessary supplies in one package
 - Step 1 and step 2 are both contained in a single outer wrapping
 - Promotes following best practice without having to remember to gather each supply individually
 - Aids in keeping *dirty* supplies separate from clean, sterile supplies
- ⇒ Available with BioPatch or Tegaderm CHG



To order product:

- Review financial impact with your Clinical Resource Director (supply chain representative)
- Determine the best plan for your facility and which products you will use
- Work with the Clinical Resource Director to order and transition to the new kit
- Plan communication and education for the nursing team (education materials are in progress)

REORDER NO.
2585

LSL
STERILE* Stérile
Central Line Kit w/ Biopatch

Not Made with Natural Rubber Latex

CONTENTS:

1 ea Step # 1 Label (Outside of first CSR Wrap)

1 ea Hand Sanitizer (in 2nd fold of CSR Wrap)
2 ea Face Mask
1 pr Vinyl Walleled Gloves
1 ea LSL Tape
1 ea Drape
3 ea Alcohol Prep Pad

1 ea Step # 2 Label (Outside of second CSR Wrap)

1 ea Hand Sanitizer (in 2nd fold of CSR Wrap)
1 pr Vinyl Walleled Gloves
1 ea Chloraprep® 3ml
1 ea Biopatch 10150
1 ea Tegaderm 1685
1 ea Foam Tray

Contents placed in order from bottom to top. Hand sanitizer in 2nd folds.

Warning: Use with care in premature infants or infants under 2 months of age. These products may cause irritation or chemical burns.

*Contents of inner wrap are sterile unless wraps has been opened or damaged. Single use only. Disposable. Contents of Foil Packets and Applicator have not been terminally sterilized by LSL Healthcare, Inc. Store between 20-25°C. Avoid freezing and excessive heat above 40°C.

REV: 7/24/2017

LSL Healthcare, Inc., Chicago, Illinois

REORDER NO.
2557

LSL
STERILE* Stérile
Central Line Kit w/ 1657 CHG

Not Made with Natural Rubber Latex

CONTENTS:

1 ea Step # 1 Label (Outside of first CSR Wrap)

1 ea Hand Sanitizer (in 2nd fold of CSR Wrap)
2 ea Face Mask
1 pr Vinyl Walleled Gloves
1 ea LSL Tape
1 ea Drape
3 ea Alcohol Prep Pad

1 ea Step # 2 Label (Outside of second CSR Wrap)

1 ea Hand Sanitizer (in 2nd fold of CSR Wrap)
1 pr Vinyl Walleled Gloves
1 ea Chloraprep® 3ml
1 ea Tegaderm 1657
1 ea Foam Tray

Contents placed in order from bottom to top. Hand sanitizer in last folds.

Warning: Use with care in premature infants or infants under 2 months of age. These products may cause irritation or chemical burns.

*Contents of inner wrap are sterile unless wraps has been opened or damaged. Single use only. Disposable. Contents of Foil Packets and Applicator have not been terminally sterilized by LSL Healthcare, Inc. Store between 20-25°C. Avoid freezing and excessive heat above 40°C.

Rev: 6/20/2017

LSL Healthcare, Inc., Chicago, Illinois

Procedure Guide for the Two Step Kit Central Line Maintenance



- Obtain the individually wrapped kit with Step 1 and Step 2 individually wrapped inside.
- Prepare the patient and provide education about the procedure.

Step 1

1. Open sterile kit (using the outer towel as a work field)
2. Sanitize hands with the individual packet of hand sanitizer in the kit
3. Don Mask
4. Mask and drape patient
5. Don sterile exam gloves
6. Stabilize the catheter, remove the old StatLock, and remove the dressing by pulling towards the insertion site (tape is provided to stabilize the line during this process and alcohol pads for cleaning the skin during removal).
7. Remove dirty gloves and prepare to use Step 2 kit



Step 2 (open cautiously as your hand sanitizer is in the fold of the kit)

1. Open sterile kit (use the outer towel as a sterile field)
2. Sanitize hands and drop StatLock onto the sterile field (obtained separate from the kit)
3. Don sterile exam gloves
4. Frictional CHG scrub of insertion and statlock site (30 second scrub and 30 second dry)
5. Apply new StatLock
6. Apply Biopatch as indicated by facility protocol
7. Apply new transparent dressing (Tegaderm)
8. Date and time the dressing (The Tegaderm has a peel away label within the border for dating and timing the dressing. You can mark on this label with a standard pen and replace on the dressing after you have completed the sterile process and washed your hands).



Note: The kit does not include a StatLock. This item will need to be obtained as a separate item.