



Hold, Transport, Store, Refrigerate, Incubate, and Dispose of Test Tubes, All in One Convenient Unit!

**test n' toss**  
Disposable Test Tube Rack



## Low-Cost Convenience:

### Spill-Resistant and Stacks for Storage, Incubator, or Refrigerator

The Test n' Toss™ is perfect for labs, hospitals, EMS, tissue recovery and organ procurement teams, physician offices and alternate care sites.

The Test n' Toss holds up to (80) 16 x 100 mm tubes in a leak-resistant plastic shell to contain spills, and is inexpensive enough to dispose of when contaminated. Clear walls allow you to see rack contents.

#### Features & Benefits:

- Economical
- Disposable
- Leak-Resistant
- Transportable
- See-Through
- Refrigeratable
- Closable
- Incubatable
- Stackable



INSERTS: 100% RECYCLED FIBERS 35% MINIMUM POST CONSUMER CONTENT. RACKS: MADE FROM RECYCLED POLYESTER.

MADE IN THE  
**USA**

## Transport / EMS:



Use a second shell as a cover to completely enclose contents for transport or storage, allowing racks to be stacked. Divider shelves can be used singly or in pairs for extra tube stability.

#### Test N' Toss™ Disposable Test Tube Rack:

Order No.:	18710
Dimensions:	
Open Rack:	11" L x 9" W x 3" H
Covered Rack:	11" L x 9" W x 5" H
Shipping Weight:	3 Lbs.
Quantity:	Carton of 10 half-shells and dividers, makes 10 open racks or 5 covered, stackable racks
Description:	Each rack holds (80) 16 x 100 mm test tubes

You'll wonder...  
how you ever managed *your lab* without us.<sup>SM</sup>

