

**SUGGESTED FORMULA**

Morphine Sulfate 25mg Slow Release Rectal Suppositories

Version number: 1.0

Volume or quantity: #100 Suppositories

Morphine Sulfate, USP (M1167)	2.5gm
Alginic Acid, NF (AL172)	25gm
Supposi-Plex R36 Suppository Base (550-83350)	Q.S.

SUGGESTED COMPOUNDING PROCEDURES

*Prior to beginning, it is necessary that any mold used be calibrated with the formula to determine the quantity of base required.

1. Calculate the required quantity of each ingredient for the total amount to be prepared
2. Accurately weigh and/or measure each ingredient
3. Pass the Alginic Acid through a #200 Mesh sieve
4. Using an appropriate heat source, melt the required quantity of Supposi-Plex R36
5. Sprinkle the Alginic Acid on the melt and mix well
6. Slowly sprinkle the Morphine Sulfate on the mixture and stir until mixed
7. Place the mixture on an ultrasonic bath for 10 minutes for thorough and complete mixing
8. Cool slightly, then pour into the mold in a continuous motion, if appropriate for the mold being used
9. Cool, trim if necessary, package and label
10. Suggested Quality assessments
 - a. % variability & % deviation from theoretical <10% - (remake if >10%)
 - b. Surface appearance and feel
 - c. Melting test
 - d. Quantity
 - e. Label - auxiliary labels, storage, BUD, compounded medication

Package in air tight, light resistant containers**Store at room temperature**

No claims are made as to the safety or efficacy of this preparation. This formulation is provided solely at the unsolicited request of the pharmacist.

Beyond-Use Dates of preparations are conservative estimates from reference books, peer-reviewed literature, and intended duration of therapy, formulation from commercially available products, organoleptic observations and current USP guidelines. Compounders may have stability studies performed by a reputable laboratory if they wish to extend the Beyond-Use Date.

It is recommended that you follow USP <795> recommendations for potency testing.

Beyond-Use Date should be assigned based on the current USP Guidelines.

Precautions:

Precautions should be taken to prevent cross-contamination and exposure of ingredients to the compounder and contamination of the preparation by the compounder. Wear appropriate protective equipment. Use safety enclosures (hoods) when weighing and mixing

Although much attention has been paid to ensure the accuracy of the formulation contained here, Spectrum Pharmacy Products accepts no liability for the loss or damage arising from reliance on the information. Compounding pharmacists using this formula take full responsibility for the formulations and hold Spectrum Pharmacy Products and Spectrum Chemical Mfg Corp. and its officers, directors and employees harmless for any claim arising from use of or reliance on information contained therein.

9/18JD

Phone: 800.370.6231 | Fax: 732.608.5420 | Internet: SpectrumRx.com | Email: sales@spectrumrx.com
© 2015 Spectrum Pharmacy Products. All rights reserved.

spectrum
PHARMACY PRODUCTS