


## N-EVAP Accessories

### Needles and Pipette Adaptors

<p>Stainless steel needles, 4in x 19 gauge, 304 grade shaft, brass plated hub</p> <p>Set of one dozen Set of 50</p>	<p><a href="#">NA0603</a> <a href="#">NA1303</a></p>	<p>\$46.00 \$182.00</p>	
<p>Stainless steel needles, 4in x 19 gauge, 304 grade shaft, plastic hub</p> <p>Set of one dozen Set of 50</p>	<p><a href="#">NA0604-19</a> <a href="#">NA1304-19</a></p>	<p>\$24.00 \$92.00</p>	
<p>FEP coated needles, 4in x 19 gauge, for corrosive solvents</p> <p>Set of one dozen Set of 50</p>	<p><a href="#">NA0603-T</a> <a href="#">NA1303-T</a></p>	<p>\$72.00 \$288.00</p>	
<p>Pasteur pipette adapter and flow controller, male luer to fluted nut</p> <p>Set of one dozen</p>	<p><a href="#">NA0636</a></p>	<p>\$93.00</p>	
<p>Pasteur pipette adapter and flow controller, male luer to 1/8 hose barb</p> <p>Set of one dozen</p>	<p><a href="#">NA0637</a></p>	<p>\$70.00</p>	

### Accessories

<p>Nitrogen Gas Generator</p> <p>20 lpm up to 48 sample positions</p>	<p>NA1955</p>	<p>\$5980.00</p>	
<p>Nitrogen Generator</p> <p>All in one nitrogen generator with compressor up to 40 lpm</p>	<p>NA1990</p>	<p>\$12,950</p>	
<p>Pressure reducing regulator for gas, 30-150 psig in, 0-30 psig out</p>	<p><a href="#">NA0630</a></p>	<p>\$112.00</p>	
<p>12 position sample holder for small glassware 5-16 mm in diameter</p> <p>Compatible with N-EVAP 111</p>	<p>NA1132</p>	<p>\$168.00</p>	
<p>6 position sample holder for large glassware 30-56 mm in diameter</p> <p>Compatible with N-EVAP 111</p>	<p>NA1113</p>	<p>\$146.00</p>	
<p>8 position sample holder for large glassware 30-75 mm in diameter</p> <p>Compatible with N-EVAP 112</p>	<p>NA1228</p>	<p>\$318.00</p>	

<p>Aluminum beads for dry bath</p> <p>3.5 liters for N-EVAP 111 5.5 liters for N-EVAP 112</p>	<p><a href="#">NA1119</a> <a href="#">NA1219</a></p>	<p>\$105.00 \$178.00</p>	
<p>Glass beads for acid resistant dry bath</p> <p>3.5 liters for N-EVAP 111-RT 5.5 liters for N-EVAP 112-RT</p>	<p><a href="#">NA1120</a> <a href="#">NA1220</a></p>	<p>\$42.00 \$73.00</p>	