



Datt Mediproducts Private Limited



VEL NeXT™

Negative Pressure Wound Therapy System



PROPERTIES

- A novel selection range of three therapies: Continuous, Intermittent and Sinusoidal
- Graphical Liquid Crystal Display (GLCD)
- Light weight device (1.5 kg approx.), easily portable
- Smart Battery Backup
- Virtual silent operation
- Well-developed Alert system with a range of Safety alarms (both audio & visual)
- Pressure Settings: A pressure setting range of 0 to -250 mmHg for three therapies
- Discrete & disposable Canister with a capacity of 500 ml and a Super absorbent polymer (SAP)
- Four different Dressing kits are available as per wound size-S, M, L, XL

INDICATIONS

- Diabetic & Neuropathic ulcer
- Acute / Chronic Wounds
- Traumatic Wounds
- Surgical Wounds
- Dehisced Wounds
- Skin Grafting
- Pressure ulcer (stage 3 & 4)
- Partial thickness burns
- Venous leg ulcer

SIZES

Each dressing kit contains a Smart Gauze, Non-adherent PU Foam, Transparent PU Film, Track Pad and a Wound measuring scale.

Small Dressing Contents:

- PU Foam 10cm x 8cm

Medium Dressing Contents:

- PU Foam 18cm x 13cm

Large Dressing Contents:

- PU Foam 26cm x 16cm

Extra Large Dressing Contents:

- PU Foam 32cm x 20cm

COLOR

- White & Light Sea Green

TEMPERATURE

- Not to exceed 25°C

SHELF LIFE

- 3 years from the date of mfg.

MATERIAL

- Foam is made of Polyester based Polyurethane (PU) Material, Transparent PU Film with adhesive

Datt Mediproducts Pvt. Ltd.
 Delivering Next-Gen Healing Solutions Globally
 www.dattmedi.com
 © All Rights Reserved

Email: advancedevice@dattmedi.com
 Phone: +91-11-47191777
 Address: 56, Community Centre,
 East Of Kailash,
 New Delhi, India - 110065

An ISO 9001 & 13485 Certified Company
 Follow Us and Stay Tuned
 @ dattmediproducts

