




Datt Mediproducts Private Limited



VELSEAL™ - PATCH

The Optimal Haemostatic Patch for Post Renal Dialysis

Patent No. US 10,124,084 B2



PROPERTIES

- A biocompatible polymeric haemostatic dressing that aids in blood clotting at the dialysis puncture site
- Highly porous
- Absorbs fluid 12-15 times of its weight
- Prevents the loss of natural clotting factors to initiate the clotting cascade
- Adhesive is non-irritant and gentle to the skin
- Biocompatible dressing
- No secondary dressing is required
- Even effective in patients with high blood pressure and/or are on anti-coagulants
- Used for haemodialysis patients

INDICATIONS

- Intended to be used as a haemostatic dressing for Post Renal Dialysis patients.

SIZES

- 5cm x 7cm

COLOR

- White

TEMPERATURE

- 8°C - 10°C

SHELF LIFE

- 2 years from the date of mfg.

MATERIAL

- Polymeric haemostatic pad on a fabric coated with adhesive

