




Datt Mediproducts Private Limited



VELSEAL™ -T

The Action Dressing

Patent No. US 10,124,084 B2

PROPERTIES

- A novel absorbent dressing containing clot promoting & clot stabilizing agents
- Used on external wounds to stop even severe arterial or venous bleeding within 4-6 minutes
- Highly porous
- No secondary bleeding
- Does not cause exothermic reactions
- Non-cytotoxic, non-irritant & non-sensitizer
- Biocompatible
- Simple to remove
- Stops bleeding even on patients with Coagulopathy and bleeding disorders (Haemophilia)
- No excretion of active compounds
- Less chances of Disseminated Intravascular Coagulation (DIC)

INDICATIONS

- Post Angiography/Angioplasty
- Central Venous Catheter insertion
- Pacemaker insertion
- Skin grafting
- Femoral arterial punctures
- Accidents, Trauma and emergency
- Amputation instead of cauterization
- Maxillo-facial trauma
- Haemophilia (External wounds only)
- Head and Neck surgeries
- Dental surgeries

SIZES

- 3cm x 3cm
- 3cm x 5cm
- 5cm x 5cm
- 7cm x 7cm
- 10cm x 10cm

COLOR

- White

TEMPERATURE

- 2°C - 8°C (Do not freeze)

SHELF LIFE

- 2 years from the date of mfg.

MATERIAL

- A Novel Polyelectrolyte complex (PEC) of biopolymers impregnated with clot promoting & stabilizing agents

