



Datt Mediproducs Private Limited

(An ISO 13485 Certified Company)

Delivering Next-gen Healing Solution Globally



DATT MEDIPRODUCTS PVT. LTD.

Datt Mediproducts Pvt. Ltd. was established in 1995, as a part of a well diversified and sound conglomerate with wide ranging business interests in India and overseas. It was initially set up as a 100% Export Oriented Unit (EOU), which commenced commercial production of a varied range of textile based medical disposables in early 1997. Now one of the largest modern plants of its kind in India, it is situated 60 km away from New Delhi in Sohna, near to Gurgaon.

In order to achieve the demanding quality standards expected in the International markets all products are made and packed in a controlled environment, which has been designed as per Class 10,000 & 100,000 specifications. Technical experts oversee the critical areas of production and packing to ensure an uncompromising adherence to quality control.

With a government recognised (DSIR) research & development lab, we committed to constant research and development to evolve new products in advanced wound management. Catering to the ever-changing needs of the world's most discriminating customers we always strive to enhance product quality.

DMPL is an ISO 13485 certified, GMP Approved unit and products are CE certified. Exporting our products to major companies in countries like Europe, Australia, New Zealand, Middle & Far-East and South Africa, we are successful in building a strong foundation globally.





VELSOFT®

Short Stretch Compression Bandage

Properties	<ul style="list-style-type: none">• Type 3(a) - Light compression bandage• 100% Cotton threads, very soft and comfortable to the skin• No rubber threads generates comfortable and safe pressure• Sealed edges, no loose threads• Latex Free• Non slippery texture• With clips• Washable, can be sterilized by steam/ ETO if required
------------	--

Colour	Flesh
--------	-------

Material	100% Cotton
----------	-------------

Indication	<ul style="list-style-type: none">• Support bandage in sprains & strains• In lymphatic disorders• Securing bandage for skin grafts• Pressure dressing for burns
------------	--

Sizes	5cm x 4m 8cm x 4m 10cm x 4m 15cm x 4m
-------	--





VELSOFT[®]-V

Properties	Elastic Crepe Bandage with Velcro Closure <ul style="list-style-type: none">• Ideal for light compression• Easy to use with Velcro closure at both ends• No clips or tapes required to hold the bandage in place• Very soft, comfortable and breathable to skin• Cotton/poly blend yarn offers excellent elasticity and fit• Latex free• Sterile (ETO/ Gamma) & Non-sterile• Individually wrapped; Sterile bandages are provided in a peel-pouch
Colour	Flesh
Material	Nylon, Viscose & Nylon coated lycra
Indication	<ul style="list-style-type: none">• Compression bandage in strains/ sprains, lymphedema & varicose vein cases• Securing bandage for skin grafts• Pressure dressings• In OT, as compression bandage
Sizes	5 cm x 4.5 m 8 cm x 4.5 m 10 cm x 4.5 m 15 cm x 4.5 m (Available in Jar pack (Non-sterile) & Peel off Pouch pack (Sterile))





VELSOFT®-BLUE

Short Stretch Compression Bandage	
Properties	<ul style="list-style-type: none">• Type 3(a) - Light compression bandage• 100% Cotton threads, very soft and comfortable to the skin• No rubber threads generates comfortable and safer pressure• Sealed edges, no loose threads• Latex Free• Non slippery texture• With clips• Washable, can be sterilized by steam/ ETO if required
Colour	Blue
Material	100% Cotton
Indication	<ul style="list-style-type: none">• Support bandage in sprains & strains• In lymphatic disorders• Securing bandage for skin grafts• Pressure dressing for burns
Sizes	5cm x 4m 8cm x 4m 10cm x 4m 15cm x 4m





VELKOMP[®]

Short Stretch Bandage for Moderate Compression	
Properties	<ul style="list-style-type: none">• Type 3(b)- Moderate compression bandage• 100% Cotton threads, very soft and comfortable to skin• No rubber threads, generates comfortable and safer pressure• Sealed edges, no loose threads• Latex free• Non slippery texture• With Clips• Washable, can be sterilized by steam/ ETO if required
Colour	Flesh Colour
Material	100% Cotton
Indication	<ul style="list-style-type: none">• Treatment of Varicose Veins, CVI, DVT• Support bandage in Sprain & Strain• Retention of dressing, splints, hot/cold packs• Prophylaxis against sport injuries
Sizes	5 cm x 4.5 m 8 cm x 4.5 m 10 cm x 4.5 m 15 cm x 4.5 m





VELFORMTM

Woven Conforming Cohesive Bandage

Properties	<ul style="list-style-type: none">• 100% Cotton, light stretch & regain• Latex free coating for cohesive effect• Bandage sticks to itself and not to the skin• Cohesiveness prevents slippage problem
Colour	Pure White
Material	100% Cotton, Latex Free Adhesive
Indication	Support underlying dressing at high mobility joints
Sizes	5 cm x 3.5 m 7.5 cm x 3.5 m 10 cm x 3.5 m 15 cm x 3.5 m





VELFORM™-C

Woven Conforming Retention Bandage	
Properties	<ul style="list-style-type: none">• Woven conforming retention bandage• 100% cotton, prevents slippage problem• Bandage can be applied on most contoured body parts• Very smooth and comfortable, no skin irritation
Colour	Pure White
Material	100% Cotton
Indication	Support underlying dressing at high mobility joints
Sizes	5 cm x 3.5 m 7.5 cm x 3.5 m 10 cm x 3.5 m 15 cm x 3.5 m





M-CREPE

Cotton Crepe Bandage	
Properties	<ul style="list-style-type: none">• Short stretch 100% Cotton crepe bandage• Type 2 bandage• Good elasticity• Allows normal skin to breathe• Soft to feel• Washing in warm water restores elasticity, bandage can be reused• With Clips
Colour	Flesh
Material	100% Cotton
Indication	<ul style="list-style-type: none">• Prevents edema formation• Support bandage in mild sprains & strains
Sizes	5 cm x 4 m 7.5 cm x 4 m 10 cm x 4 m 15 cm x 4 m





VELSOFT[®]-PLUS

Woven Cohesive Moderate Compression Bandage	
Properties	<ul style="list-style-type: none">• 100% Cotton woven bandage• Coated with latex free components• Provides light to moderate compression with Cohesive effect- sticks to itself• Flexible substrate allows it to mould to most difficult of body contours• Hand tearable in both directions- easy to use• Air permeable- open construction allows skin to breathe
Colour	Blue
Material	Cotton, Acrylic Adhesive, Viscose
Indication	<ul style="list-style-type: none">• Retention of dressing, splints, hot/cold packs• Prophylaxis against sports injuries• Useful in sprains & strains
Sizes	10 cm x 2 m (Unstretched length)





COKOMPTM

Non Woven Cohesive Moderate Compression Bandage

Properties	<ul style="list-style-type: none">• Cohesive on both sides, coated with latex free components, adheres only to itself, not to skin• Hand tearable material• Non-slippery in texture• Provides sustained & reliable compression upto 23mmHG
Colour	Brown
Material	Non Woven, Elastic fibers, Latex free acrylic adhesive
Indication	<ul style="list-style-type: none">• Treatment of Sprains & Strains, CVI, DVT, Varicose veins & Venous leg ulcers• Overwrap for splints• For securing dressings
Sizes	5 cm x 4.5 m 7.5 cm x 4.5 m 10 cm x 4.5 m 15 cm x 4.5 m





VELSOFT[®]-COHESIVE

Short Stretch Cohesive Compression Bandage	
Properties	<ul style="list-style-type: none">• Type 3(b)- Moderate compression bandage• 100% Cotton threads, very soft and comfortable to skin• No rubber threads, generates comfortable and safe pressure• Sealed edges, no loose threads• No slippage problem• No clips are required• Bandage sticks to itself and not to the skin
Colour	Flesh
Material	100% Cotton
Indication	<ul style="list-style-type: none">• Useful in Varicose Veins, CVI, DVT• Support bandage in Sprains & Strains• As securing bandage for skin grafts• Pressure dressing in burns• In lymphatic disorder & to eliminate venous statis
Sizes	5 cm x 4.5 m 8 cm x 4.5 m 10 cm x 4.5 m 15 cm x 4.5 m





VELSOFT[®]-ADHESIVE

Elastic Adhesive Bandage	
Properties	<ul style="list-style-type: none">• 100% Latex free, no Latex allergy & Skin irritation• Leaves no traces of adhesive on the skin after removal• Skin friendly• Excellent adhesion• Porous adhesive, permeable to air & moisture• Prevents skin maceration• Offers superior stretch
Colour	Flesh, White adhesive
Material	100% Cotton, Zinc Oxide
Indication	<ul style="list-style-type: none">• Post operative dressing• For compression & support of muscles and joint• Angioplasties
Sizes	8 cm x 1 m 10 cm x 1 m 10 cm x 2/3 m 6 cm x 4/6 m 8 cm x 4/6 m 10 cm x 4/6 m 15 cm x 4/6 m





VELCOOL[®]

Combine Cooling & Compression

Properties	<ul style="list-style-type: none">• Perfectly adheres to R.I.C.E Method of treatment• Cold bandage combines cooling and compression for treating pains, aches, sprains & strains• Provides pain relief in seconds and cools for hours• Bandage contains a specially formulated solution that begins to evaporate upon contact with skin thereby cooling the applied area• Keeps the applied area cool for 4-6 hours• Quickly reduces swelling & pain• After 6 hours the bandage works as a Compression Bandage• Cohesive bandage retains the cold bandage and provides complete and comprehensive cooling effect
Colour	Cooling bandage- Flesh Compression bandage- Blue
Material	Cooling bandage & Compression bandage is made of 100% cotton, Latex free acrylic adhesive coating on cohesive bandage
Indication	Most beneficial treatment of sprains & strains in the first 24 Hrs - 72 Hrs after sustaining injury
Sizes	Impregnated Bandage 10 cm x 5 m Cohesive Bandage 10 cm x 2 m





VELSOFT[®]-FR

Fixing Roll for Wide Area Fixation of Wound Dressing

Properties

- Non-Woven Fixing Roll
- 100% Latex Free
- Skin Friendly
- Hypoallergenic adhesive
- Self adhesive
- Permeable to air and water vapour
- Radio transparent
- Allows joints to be dressed without constraining movement
- Adapts well even on high contoured parts of the body

Colour

White

Material

Non Woven, Acrylic Adhesive

Indication

- Wide area fixation of wound dressing
- Secure fixation of tubes & catheters

Sizes

5 cm x 2 m
10 cm x 2 m
15 cm x 2 m
5 cm x 10 m
10 cm x 10 m
15 cm x 10 m





VELCAST™

Properties	Orthopaedic Fiberglass Casting Tape <ul style="list-style-type: none"> • Excellent strength, rigidity & durable • Good setting time & smooth moulding • Excellent X-ray radiolucency • Easy to apply and remove • Lighter and stronger than the traditional POP bandage • Water resistant • Conforms to the area of body smoothly without wrinkles • High quality knitted substrate and special resin protects against unexpected sudden impact • Functional strength achieved in 20 minutes. Cost effective
Colour	White, Pink, Blue & Green (5 cm cast is available in polyester material with graphics)
Material	Fiberglass, Polyester
Indication	For Immobilization
Sizes	5 cm x 3.6 m 7.5 cm x 3.6 m 10 cm x 3.6 m 12.5 cm x 3.6 m 15 cm x 3.6 m





VELPOP®

Properties	Plaster of Paris Bandage
	<ul style="list-style-type: none">• Combines fast setting, high quality Plaster of Paris with conformable and stable interlocked weave of the best quality cloth• The rich creamy texture allows easy handling• Interlocked weave of the cloth helps in plaster retention & provides high degree of accuracy• High quality creamy textured Plaster of Paris makes moulding easy• Moisture resistant packaging improves the protection of bandages prior to use• Setting time is just 2-4 minutes• Excellent whiteness• Safe for skin
Colour	White
Material	Woven Cloth, Plaster of Paris
Indication	For immobilization to protect limbs with broken bones
Sizes	5cm x 2.7 m 7.5cm x 2.7 m 10cm x 2.7 m 15cm x 2.7 m 20cm x 2.7 m





VELNETTM

Properties	100% Cotton Stockinette
	<ul style="list-style-type: none">• Protects the skin surface under compression• 100% cotton for non-irritative comfort• Easy application• Light weight and strong protection• Conforms to joints, will not wrinkle or slip• Durable and economical, may be washed for reuse• Extra protection for frequent bandage changes• Highly absorbent and well tolerated• Steam autoclavable for sterile applications
Colour	White
Material	100% cotton
Indication	Protect skin from friction under plaster/fibreglass bandage
Sizes	5cm x 1m 7.5cm x 1.5m 10cm x 1.5m 5cm x 10m 7.5cm x 10m 10cm x 10m





VELROLIL[®]

Cast Padding

Properties	<ul style="list-style-type: none">• Synthetic cast padding provides padding around the bony structure• Provides comfort and protection• Uniform distribution of pressure• Air permeable• Protects sensitive skin
------------	--

Colour	Pure White
--------	------------

Material	Viscose
----------	---------

Indication	Padding around bony structure
------------	-------------------------------

Sizes	10 cm x 3 m 15 cm x 3 m
-------	----------------------------





SOFTPACK[®]

	Sterile Dressing Kits
Properties	<ul style="list-style-type: none">• Ready to use sterile dressing kits• Four different kits are available with different contents• Sterilized by Gamma Irradiation/ ETO
Indication	Cleaning of wounds
Content	<p>Kit 1 : (Small)</p> <ul style="list-style-type: none">• Dressing forceps - 1• Tray with Gallipots - 1• Cotton Balls- 2• Gauze swab (3" x 3" x 12 ply) - 3• Absorbent dressing pad (3" x 7.5") - 1• Hand wipe tissue - 1 <p>Kit 2 : (Large)</p> <ul style="list-style-type: none">• Sterile Drape - 1• Dressing forceps - 1• Tray with Gallipots - 1• Cotton Balls- 4• Gauze swab (3" x 3" x 12 ply) - 6• Absorbent dressing pad (3" x 7.5") - 1• Hand wipe tissue - 1 <p>Kit 3 :</p> <ul style="list-style-type: none">• Cotton Balls- 4• Gauze swab (4" x 4" x 12 ply) - 4• Absorbent dressing pad (3" x 7.5") - 1• Hand wipe tissue - 1• Sterile sheet - 1 <p>Kit 4 :</p> <ul style="list-style-type: none">• Dressing forceps - 1• Tray with Gallipots - 1• Cotton Balls- 2• Gauze swab (3" x 3" x 8 ply) - 3• Absorbent dressing pad (3" x 7.5") - 1• Examination Gloves - 1 pair• Hand wipe tissue - 1• Sterile sheet - 1
Also Available	<ol style="list-style-type: none">1. Catheterization Kit2. Dialysis Kit3. Customer Specific Kit





SOFTSWAB[®]

Properties	Sterile Gauze Swabs E.P <ul style="list-style-type: none">• Made of 100% Cotton Gauze Fabric• Conform to EP- 17 Standards• All swabs are with folded edges (American Fold), therefore no loose threads• Available with or without X-ray detectable thread• Bleaching done by using Hydrogen Peroxide and de-mineralized water• ETO/ Gamma sterilization
Colour	White
Material	100% Cotton
Indication	<ul style="list-style-type: none">• In OT to absorb blood• Wound dressing
Sizes	5cm x 5cm x 8/12/16 ply 7.5cm x 7.5cm x 8/12/16 ply 10cm x 10cm x 8/12/16 ply





VELSWAB[®]

Properties	Unsterile Gauze Swabs E.P <ul style="list-style-type: none">• Made of 100% Cotton Gauze Fabric• Conforms to E.P-17 Standards• All Swabs are with folded edges (American Fold), therefore no loose threads• Available with or without X-ray detectable thread• Bleaching done by using Hydrogen peroxide and de-mineralized water
Colour	White, Green
Material	100% Cotton
Indication	<ul style="list-style-type: none">• In OT to absorb blood• Wound dressing
Sizes	5cm x 5cm x 8/12/16ply 7.5cm x 7.5cm x 8/12/16ply 10cm x 10cm x 8/12/16ply Pack of 100's or Steam compatible pack of 2's/4's/5's with or without sterilization indicator or as desired by the Customer





SOFTLAPS[®]

Properties	Laparotomy Sponge <ul style="list-style-type: none">• Made of 100% Cotton Gauze Fabric• Bleaching done by using Hydrogen Peroxide & de-mineralized water• All four edges are sealed to prevent lint residue• Available with or without X-ray detectable thread• Tape stitched at the corner• Available in sterile and non-sterile packaging• Tied in bundles of 5's with red thread for easy counting
Colour	White
Material	100% Cotton
Indication	In OT to absorb Blood
Sizes	20 cm x 20 cm x 4/6/8/12/16 ply 25 cm x 25 cm x 4/6/8/12/16 ply 30 cm x 30 cm x 4/6/8/12/16 ply 25 cm x 40 cm x 4/6/8/12/16 ply 40 cm x 40 cm x 4/6/8/12/16 ply 45 cm x 45 cm x 4/6/8/12/16 ply





SOFTPADTM

	Sterile Combine Dressing Pad
Properties	<ul style="list-style-type: none">• Combine dressing of cotton pad with absorbing overwrap of non-woven material• Individually packed in temper proof wraps• Sterilized by Gamma irradiation/ ETO• Non-adherent to the wound• Highly absorbent cotton pad ensures slow strike-through even in case of heavy exudates
Colour	Pure White
Material	Non Woven, Polyester
Indication	Wound dressing used to manage heavy exudates
Sizes	7.5cm x 7.5cm 10cm x 10cm 10cm x 15cm 10cm x 20cm





SOFTPADTM

Eye Pad with Adhesive Borders

Properties	<ul style="list-style-type: none">• Non-woven, Latex Free adhesive eye pad with adhesive borders• Acrylic based Hypoallergenic adhesive• Unique non-adherent absorbent pad coated with polyethylene mesh won't stick to the wound• Gamma sterilization
Colour	White
Material	Non Woven, Acrylic adhesive
Indication	Eye injury care
Sizes	10 cm x 7.2 cm





VELPORE™

Properties	Hypoallergenic & Latex free Surgical Paper Tape
	<ul style="list-style-type: none">• Hypoallergenic and 100% Latex free• White, lightweight, non-woven paper tape• Porous, allowing skin to breath• Gentle to skin• Holds well on skin• Conformable• Water repellent• Easily hand tearable• No traces of adhesive remains on the skin after removal
Colour	White
Material	Non Woven, Acrylic Adhesive
Indication	Perfect for securing of: - Dressing - Tubes - Catheters, etc.
Sizes	1/2" x 9.14/5 m 1" x 9.14/5 m 2" x 9.14/5 m 3" x 9.14/5 m Also available with Dispenser pack in 9.14 m length Available in Individual box packing by the brand name of Velpore- SR





VELPORE-T™

	Transparent Perforated Surgical PE Tape
Properties	<ul style="list-style-type: none">• Micro-porous structure ensures high air & vapour permeability• Prevents skin from maceration• Thin & excellent adhesive- comfortable for the patients• Easy to handle with gloves
Colour	White
Material	Transparent, acrylic adhesive
Indication	<ul style="list-style-type: none">• Fixation of dressings on sensitive skin• Securing of : - I.V Tubing<ul style="list-style-type: none">- Dressings- Catheters, etc.
Sizes	1/2" x 9.14 m/ 5 m 1" x 9.14 m/ 5 m 2" x 9.14 m/ 5 m 3" x 9.14 m/ 5 m





VELSILK™

Properties	Hypoallergenic & Latex free Silk-like Tape <ul style="list-style-type: none">• Latex-free, hypoallergenic silk-like tape• Adheres well with good initial stick, High strength• Water resistant
Colour	White
Material	Silk-like cloth, Acrylic adhesive
Indication	<ul style="list-style-type: none">• Fixation of dressings on sensitive skin• Securing of : - I.V Tubing<ul style="list-style-type: none">- Catheters- Probes- Cannulae• Stabilizing finger splints
Sizes	1/2" x 9.14 m 1" x 9.14 m 2" x 9.14 m 3" x 9.14 m





VELFIX[®]-T

Properties	Transparent Sterile I.V Dressing <ul style="list-style-type: none">• Offers a transparent, waterproof film that serves as a bacterial barrier yet is breathable• 30 µ thick PU film ensures increased breathability & softness• Unique coat technology with patterned adhesive offers high MVTR of 2685 GSM/24 Hrs• Latex free & hypoallergenic adhesive that gently, yet securely adheres to skin• Comes with 2 securing tapes & 1 paper strip to mention date of application• Gamma Sterile
Colour	Transparent
Material	Polyurethane film, Acrylic adhesive
Indication	Sterile I.V dressing in cannula fixation
Sizes	6 cm x 5 cm 6 cm x 7 cm 6 cm x 8.5 cm 7 cm x 9 cm 8.5 cm x 10 cm 10 cm x 12 cm





VELFIX[®]-T+PAD

Properties	Transparent film dressing with Non-adherent pad <ul style="list-style-type: none">• Offers a transparent, waterproof film that serves as a bacterial barrier yet is breathable• 30 µ thick PU film ensures increased breathability & softness• Unique coat technology with patterned adhesive offers high MVTR of 2685 GSM/24 Hrs• Latex free, non-adherent pad does not adhere to the wound bed, reducing pain upon removal• Easy to mould around corners and is patient & clinician friendly• Gamma Sterile
Colour	Transparent
Material	Polyurethane film, Absorbent Pad, Acrylic adhesive
Indication	<ul style="list-style-type: none">• Post operative wounds• Acute wounds (Abrasions & Laceration)• Small or medium exudate level chronic wounds
Sizes	5 cm x 7.5 cm 6 cm x 10 cm 9 cm x 12.5 cm 8 cm x 15 cm 10 cm x 10 cm 10 cm x 20 cm 10 cm x 25 cm 10 cm x 30 cm 10 cm x 35 cm





VELFIX[®]-EASY



Properties	Non-Woven sterile I.V Dressing <ul style="list-style-type: none"> • Sterile, Non-woven, IV dressing comes with two securing tapes • Latex Free hypoallergenic adhesive used • Highly porous & air permeable • Protects injection site against dust, germs and other outside contaminants
Colour	White
Material	Non-woven, acrylic adhesive
Indication	I.V cannula fixation
Sizes	6 cm x 7 cm (S) 7 cm x 9 cm (M) 8.5 cm x 10 cm (L)





VELFIX®-EASY+PAD

Properties	Non-woven dressing with Non-adherent pad <ul style="list-style-type: none">• Unique non-adherent absorbent pad coated with polyethylene mesh won't stick to the wound and absorbs small amounts of exudate, allowing normal wound healing process to occur with less pain and trauma• Micropores in the non-woven material structure provides an adequate rate of moisture vapour permeability• Covered with hypoallergenic acrylic adhesive• ETO Sterile
Colour	White
Material	Non-woven, acrylic adhesive
Indication	<ul style="list-style-type: none">• Post operative wounds• Acute wounds (Abrasions & Laceration)• Small or medium exudates level chronic wounds
Sizes	5 cm x 7.5 cm 6 cm x 10 cm 10 cm x 10 cm 9 cm x 12.5 cm 8 cm x 15 cm 10 cm x 20 cm 10 cm x 25 cm 10 cm x 30 cm 10 cm x 35 cm





SILVEL™

Nanocrystalline Silver coated Antimicrobial Barrier Dressing	
Properties	<ul style="list-style-type: none">• Three layered dressing. Wound contact layer is made of 100% Cotton, middle layer is made of Viscose and outer layer is made of polurethane net• Sustained release of Ag⁰ is 294 ppm/ 24 hours• Ag⁰ content in Silvel is 438 mg/100 cm²• Provides long lasting anti-microbial barrier (approx. 5 days)• Effective against resistance bacteria like MRSA/VRE/Pseudomonas• Maintains moist wound enviroment optimal for fast healing• Silvel needs to be moistened only once prior to use, repeated moistening is not required• Non adherent in nature• Gamma Sterile
Colour	Amber
Material	Nanocrystalline Silver coated non adherent cotton fabric, Absorbent Pad
Indication	<ul style="list-style-type: none">• Partial thickness burns• Pressure ulcer• Diabetic foot ulcer• Infected wounds• Surgical wounds• Venous leg ulcer
Sizes	5 cm x5 cm 10 cm x 10 cm 10 cm x 20 cm 15 cm x 15 cm 20 cm x 20 cm 20 cm x 40 cm 40 cm x 40 cm 10 cm x 120 cm





VELFOUR[®]

Four layer Bandaging System for Treatment of Venous Leg Ulcers & Lymphedema

Properties

- 4 layer bandaging system for treating venous leg ulcers & associated conditions
- Four kits are available which contains number of different components
- Kit A & Kit A Plus are for those patients having ankle circumference between 18-25cms, whereas Kit B and Kit B Plus are for those who have ankle circumference more than 25cms
- Provides approximately 40mmHG of High Sustained Compression Therapy (HSCT) pressure at the ankle, reducing to 17mmHg at knee
- Maintains effective level of compression even up to a full week after application
- Orthopaedic wool protects the bony structures
- Sufficient absorption capacity to manage exudates for up to a full week without needing re-application

Indication

Venous leg ulcers, Lymphedema

Content

Kit A : (18 cm- 25 cm)

- Smart Gauze: 10 cm x 20 cm
- Velfour #1: Orthopaedic wool 10 cm x 3.5 m
- Velfour #2: Cotton Crepe bandage 10 cm x 4.5 m
- Velfour #3: Light weight class 3a Elastic Compression bandage 10 cm x 8.7 m
- Velfour #4: Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)

Kit B : (>25 cm)

- Smart Gauze: 10 cm x 20 cm
- Velfour #1: Orthopaedic wool 10 cm x 3.5 m
- Velfour #2: Cotton Crepe bandage 10 cm x 7 m
- Velfour Plus: Class 3c High Compression bandage 10 cm x 5 m
- Velfour #4: Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)

Kit A Plus : (18 cm- 25 cm)

- Silvel Ultra: 10 cm x 20 cm
- Velfour #1: Orthopaedic wool 10 cm x 3.5 m
- Velfour #2: Cotton Crepe bandage 10 cm x 4.5 m
- Velfour #3: Light weight class 3a Elastic Compression bandage 10 cm x 8.7 m
- Velfour #4: Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)

Kit B Plus : (>25cm)

- Silvel Ultra: 10 cm x 20 cm
- Velfour #1: Orthopaedic wool 10 cm x 3.5 m
- Velfour #2: Cotton Crepe bandage 10 cm x 7 m
- Velfour Plus: Class 3c High Compression bandage 10 cm x 5 m
- Velfour #4: Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)



Silvel Ultra is Nanocrystalline Silver treated non-adherent antimicrobial wound contact layer. It kills the broad spectrum bacteria like MRSA/ VRE/ Pseudomonas. No need to use any antimicrobial solution, prior to use Kit A Plus & Kit B Plus





VIELNEXT™

	Negative Pressure Wound Therapy System
Properties	<ul style="list-style-type: none">• Large display screen, weight of device is just 1.4 kg (approx)• Over 18-hours (approx) battery run time• Virtual silent operation, audible & visual alarms for easy troubleshooting• Wound contact layer is made of Nano-crystalline silver• Discreet & disposable Canister with a capacity of 500 ml and Super absorbent polymer (SAP)• Three different Dressing kits are available as per wound size
Weight	1.4 Kg (Approx)
Dimension	21.8 L x 12.0 W x 14.0 H (in cm ³)
Colour	White & Light Sea Green
Material	Foam is made of Polyester based Polyurethane (PU) Material, Transparent PU Film, Acrylic adhesive
Indication	Diabetic & Neuropathic Ulcer, Acute/Chronic wounds, Traumatic wounds, Surgical wounds, Dehiscent wounds, Skin grafting, Pressure ulcer, 1 st & 2 nd degree burns, Venous leg ulcer
Sizes	<p>Small Dressing Contents:</p> <ul style="list-style-type: none">• Non adherent gauze layer with Nano crystalline silver 10cm x 10cm• Foam 10cm x 8cm• Transparent Film PU 30cm x 40cm• Track Pad <p>Medium Dressing Contents:</p> <ul style="list-style-type: none">• Non adherent gauze layer with Nano crystalline silver 12cm x 10cm• Foam 18cm x 13cm• Transparent Film PU 30cm x 40cm• Track Pad <p>Large Dressing Contents:</p> <ul style="list-style-type: none">• Non adherent gauze layer with Nano crystalline silver 15cm x 15cm• Foam 26cm x 16cm• Transparent Film PU 30cm x 40cm• Track Pad





VIELSEAL

The Action Dressing

Properties	<ul style="list-style-type: none">• Ready to use Foam based dressing for controlling traumatic bleeding• Stops arterial and venous bleeding within 4 minutes• Biocompatible and Bioabsorbable• Does not stick to the wound• Fits any wound shape, available in different sizes• Indent for use as a topical dressing
------------	---

Colour	Brown
--------	-------

Material	Polyelectrolyte complex (PEC) of biopolymers, Clot promoting and stabilizing agents
----------	---

Indication	Trauma Cases, Gunshot wounds, Accidents, Injury or Casualty
------------	---

Storage	24°C
---------	------

Sizes	7 cm x 7 cm 10 cm x 10cm
-------	-----------------------------





VIELSEAL-T

The Action Dressing	
Properties	<ul style="list-style-type: none">• Ready to use Foam based dressing• Stops arterial and venous bleeding within 4 minutes• Biocompatible and Bioabsorbable• Does not stick to the wound• Fits any wound shape, available in different sizes• Indent for use as a topical dressing
Colour	White
Material	Polyelectrolyte complex (PEC) of biopolymers, Clot promoting and stabilizing agents
Indication	After Angiography and Angioplasty, Trauma Cases, Gunshot wounds, Accidents, Injury or Casualty
Storage	8°C
Sizes	7 cm x 7 cm 10 cm x 10cm

