



Datt Mediproducts Private Limited

# COKOMP™

## Non Woven Cohesive Moderate Compression Bandage



### PROPERTIES

- Cohesive on both sides, coated with latex free component, adhere only to itself, not to skin
- Hand tearable material
- Non-slippery in texture
- Provides sustained & reliable compression upto 23mmHg

### COLOR

- Brown

### MATERIAL

- Non Woven, Elastic fibers, Latex free acrylic adhesive



### INDICATION

- Treatment of Sprains & Strains, CVI,
- DVT, Varicose veins & Venous leg ulcers
- Overwrap for splints
- For securing dressings



### SIZES

5 cm x 4.5 m

7.5 cm x 4.5 m

10 cm x 4.5 m

15 cm x 4.5 m

