



Datt Mediproducts Private Limited



DMP™ COTTON CREPE

Cotton crepe Bandage

PROPERTIES

- Short stretch 100% Cotton crepe bandage
- Type 3A bandage
- Good elasticity
- Allows normal skin to breathe
- Soft to feel
- Washing in warm water restores elasticity, bandage can be reused
- With Clips

COLOR

- Flesh

MATERIAL

- 100% Cotton



INDICATION

- Prevents edema formation
- Support bandage in mild sprains & strains



SIZES

5 cm x 4 m

7.5 cm x 4 m

10 cm x 4 m

15 cm x 4 m

Available in BOPP & Jar packing.

