



Datt Mediproducts Private Limited



# VELCOOL®

## Combine Cooling & Compression



### PROPERTIES

- Perfectly adheres to R.I.C.E Method of treatment
- Cold bandage combines cooling and compression for treating pains, aches, sprains & strains
- Provides pain relief in seconds and cools for hours
- Bandage contains a specially formulated solution that begins to evaporate upon contact with skin thereby cooling the applied area
- Keeps the applied area cool for 4-6 hours
- Quickly reduces swelling & pain
- After 6 hours the bandage works as a Compression Bandage
- Cohesive bandage retains the cold bandage and provides complete and comprehensive cooling effect

### COLOR

- Cooling bandage- Flesh
- Compression bandage- Blue

### MATERIAL

- Cooling bandage & Compression bandage is made of 100% cotton, Latex free acrylic adhesive coating on cohesive bandage



### INDICATION

- Most beneficial treatment of sprains & strains in the first 24 Hrs - 72 Hrs after sustaining injury



### SIZES

Impregnated Bandage	Cohesive Bandage
10 cm x 5 m	10 cm x 2 m

