



Datt Mediproducts Private Limited



VELFORM™-C

Woven Conforming Retention Bandage



PROPERTIES

- Woven conforming retention bandage
- 100% cotton, prevents slippage problem
- Bandage can be applied on most contoured body parts
- Very smooth and comfortable, no skin irritation

COLOR

- Pure White

MATERIAL

- 100% Cotton



INDICATION

- Support underlying dressing at high mobility joints



SIZES

5 cm x 3.5 m

7.5 cm x 3.5 m

10 cm x 3.5 m

15 cm x 3.5 m

