



Datt Mediproducts Private Limited



VELFORM™

Woven Conforming Cohesive Bandage



PROPERTIES

- 100% Cotton, light stretch & regain
- Latex free coating for cohesive effect
- Bandage sticks to itself and not to the skin
- Cohesiveness prevents slippage problem

COLOR

- Pure White

MATERIAL

- 100% Cotton, Latex Free Adhesive



INDICATION

- Support underlying dressing at high mobility joints



SIZES

5 cm x 3.5 m

7.5 cm x 3.5 m

10 cm x 3.5 m

15 cm x 3.5 m

