



Datt Mediproducts Private Limited



VELFOUR[®] PLUS

Elastic compression Bandage



PROPERTIES

- Type 3C bandage suitable for high compression requirement.
- Long stretch, soft texture, comfortable to skin.
- Offers low working pressure & high resting pressure.
- Sealed edges, no loose threads.
- Washable & can be sterilized by steam/ ETO if required.

COLOR

- Brown

MATERIAL

- Cotton, lycra



INDICATION

- Venous insufficiency, varicose veins, reduction of venous & post traumatic oedema
- As a compression bandage in strain, sprains & injuries to tendons



SIZES

5 cm x 4 m

8 cm x 4 m

10 cm x 4 m

15 cm x 4 m

