



Datt Mediproducts Private Limited



VELFOUR[®]

Four layer Bandaging System for Treatment of Venous Leg Ulcers & Lymphedema



PROPERTIES

- 4 layer bandaging system for treating venous leg ulcers & associated conditions
- Four kits are available which contains number of different components
- Kit A & Kit A Plus are for those patients having ankle circumference between 18-25cms, whereas Kit B and Kit B Plus are for those who have ankle circumference more than 25cms
- Provides approximately 40mmHG of High Sustained Compression Therapy (HSCT) pressure at the ankle, reducing to 17mmHg at knee
- Maintains effective level of compression even up to a full week after application
- Orthopaedic wool protects the bony structures
- Sufficient absorption capacity to manage exudates for up to a full week without needing re-application

COLOR

- Flesh

CONTENT

Kit A : (18 cm- 25 cm)

- **Smart Gauze:** 10 cm x 20 cm
- **Velfour #1:** Orthopaedic wool 10 cm x 3.5 m
- **Velfour #2:** Cotton Crepe bandage 10 cm x 4.5 m
- **Velfour #3:** Light weight class 3a Elastic Compression bandage 10 cm x 8.7 m
- **Velfour #4:** Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)



INDICATION

- Venous leg ulcers, Lymphedema

Kit B : (>25 cm)

- **Smart Gauze:** 10 cm x 20 cm
- **Velfour #1:** Orthopaedic wool 10 cm x 3.5 m
- **Velfour #2:** Cotton Crepe bandage 10 cm x 7 m
- **Velfour Plus:** Class 3c High Compression bandage 10 cm x 5 m
- **Velfour #4:** Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)

Kit A Plus : (18 cm- 25 cm)

- **Silvel Ultra:** 10 cm x 20 cm
- **Velfour #1:** Orthopaedic wool 10 cm x 3.5 m
- **Velfour #2:** Cotton Crepe bandage 10 cm x 4.5 m
- **Velfour #3:** Light weight class 3a Elastic Compression bandage 10 cm x 8.7 m
- **Velfour #4:** Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)

Kit B Plus : (>25cm)

- **Silvel Ultra:** 10 cm x 20 cm
- **Velfour #1:** Orthopaedic wool 10 cm x 3.5 m
- **Velfour #2:** Cotton Crepe bandage 10 cm x 7 m
- **Velfour Plus:** Class 3c High Compression bandage 10 cm x 5 m
- **Velfour #4:** Class 3b Cohesive Compression bandage 10 cm x 4.5 m. (It retains bandage position)

Silvel Ultra is Nanocrystalline Silver treated non-adherent antimicrobial wound contact layer. It kills the broad spectrum bacteria like MRSA/ VRE/ Pseudomonas. No need to use any antimicrobial solution, prior to use Kit A Plus & Kit B Plus