



Datt Mediproducts Private Limited



VELKOMP[®]

Short Stretch Bandage for Moderate Compression



PROPERTIES

- Type 3(b)- Moderate compression bandage
- 100% Cotton threads, very soft and comfortable to skin
- No rubber threads, generates comfortable and safer pressure
- Sealed edges, no loose threads
- Latex free
- Non slippery texture
- With Clips
- Washable, can be sterilized by steam/ ETO if required

COLOR

- Flesh Colour

MATERIAL

- 100% Cotton



INDICATION

- Treatment of Varicose Veins, CVI, DVT
- Support bandage in Sprain & Strain
- Retention of dressing, splints, hot/cold packs
- Prophylaxis against sport injuries



SIZES

5 cm x 4.5 m
4 cm x 5 m
10 cm x 10 m

8 cm x 4.5 m
6 cm x 5 m
12 cm x 5 m

10 cm x 4.5 m
8 cm x 5 m
12 cm x 10 m

15 cm x 4.5 m
10 cm x 5 m

