



Datt Mediproducts Private Limited



# VELSOFT® - ADHESIVE

## Elastic Adhesive Bandage



### PROPERTIES

- 100% Latex free, no Latex allergy & Skin irritation
- Leaves no traces of adhesive on the skin after removal
- Skin friendly
- Excellent adhesion
- Porous adhesive, permeable to air & moisture
- Prevents skin maceration
- Offers superior stretch

### COLOR

- Flesh, White adhesive

### MATERIAL

- 100% Cotton, Zinc Oxide



### INDICATION

- Support bandage in sprains & strains
- In lymphatic disorders
- Securing bandage for skin grafts
- Pressure dressing for burns



### SIZES

8 cm x 1 m

8 cm x 4/6 m

10 cm x 1 m

10 cm x 4/6 m

10 cm x 2/3 m

15 cm x 4/6 m

6 cm x 4/6 m

