



Datt Mediproducts Private Limited



VELSOFT® - BLUE

Short Stretch Compression Bandage



PROPERTIES

- Type 3(a) - Light compression bandage
- 100% Cotton threads, very soft and comfortable to the skin
- No rubber threads generates comfortable and safer pressure
- Sealed edges, no loose threads
- Latex Free
- Non slippery texture
- With clips
- Washable, can be sterilized by steam/ ETO if required

COLOR

- Blue

MATERIAL

- 100% Cotton



INDICATION

- Support bandage in sprains & strains
- In lymphatic disorders
- Securing bandage for skin grafts
- Pressure dressing for burns



SIZES

5 cm x 4 m

8 cm x 4 m

10 cm x 4 m

15 cm x 4 m

