



Datt Mediproducts Private Limited



# VELSOFT® - COHESIVE

## Short Stretch Cohesive Compression Bandage



### PROPERTIES

- Type 3(b)- Moderate compression bandage
- 100% Cotton threads, very soft and comfortable to skin
- No rubber threads, generates comfortable and safe pressure
- Sealed edges, no loose threads
- No slippage problem
- No clips are required
- Bandage sticks to itself and not to the skin

### COLOR

- Flesh

### MATERIAL

- 100% Cotton



### INDICATION

- Useful in Varicose Veins, CVI, DVT
- Support bandage in Sprains & Strains
- As securing bandage for skin grafts
- Pressure dressing in burns
- In lymphatic disorder & to eliminate venous stasis



### SIZES

5 cm x 4.5 m

8 cm x 4.5 m

10 cm x 4.5 m

15 cm x 4.5 m

