



Datt Mediproducts Private Limited



VELSOFT® COMPRESSION TAPE

Latex - Free



PROPERTIES

- 100% Latex free, no Latex allergy
- Acrylic adhesive minimises skin injuries
- Provides equivalent compression as compared to traditional elastic adhesive bandages.
- Excellent adhesion
- Skin friendly
- Leaves no traces of adhesive on the skin after removal
- Permeable to air & moisture
- Provides superior stretch

COLOR

- Flesh

MATERIAL

- Cotton , lycra



INDICATION

- Post operative dressing
- For compression & support of muscles and joint
- Angioplasties



SIZES

5 cm x 4/6m

10cm x 4/6m

