



Datt Mediproducts Private Limited



VELSOFT® I

Elastic Compression Bandage



PROPERTIES

- Type 3(a) - Light compression bandage
- Very soft and comfortable to the skin
- Sealed edges, no loose threads
- Latex Free
- Non slippery texture
- With clips
- Washable, can be sterilized by steam/ ETO if required

COLOR

- Flesh

MATERIAL

- Cotton , lycra



INDICATION

- Support bandage in sprains & strains
- In lymphatic disorders
- Securing bandage for skin grafts
- Pressure dressing for burns



SIZES

5 cm x 4.5 m

8 cm x 4.5 m

10 cm x 4.5 m

15cm x 4.5 m

