



Datt Mediproducts Private Limited



VELSOFT® - PLUS

Woven Cohesive Moderate Compression Bandage



PROPERTIES

- 100% Cotton woven bandage
- Coated with latex free components
- Provides light to moderate compression with Cohesive effect- sticks to itself
- Flexible substrate allows it to mould to most difficult of body contours
- Hand tearable in both directions- easy to use
- Air permeable- open construction allows skin to breathe

COLOR

- Blue

MATERIAL

- Cotton, Acrylic Adhesive, Viscose



INDICATION

- Retention of dressing, splints, hot/cold packs
- Prophylaxis against sports injuries
- Useful in sprains & strains



SIZES

10 cm x 2 m (Unstretched length)

