



Datt Mediproducts Private Limited



VELSPAN™

Strong Graduated long stretch bandage



PROPERTIES

- Type 3C bandage suitable for high compression requirement
- Long stretch, soft texture, comfortable to skin
- Offers low working pressure & high resting pressure
- Indicator marking on the bandage allows the precise control of pressure on application
- Stretch the bandage to 50% (rectangle marking becomes square)
- Sealed edges, no loose threads
- Washable & can be sterilized by steam/ ETO if required

COLOR

- Brown

MATERIAL

- Cotton, lycra



INDICATION

- Venous insufficiency, varicose veins, reduction of venous & post traumatic oedema.
- As a compression bandage in strain, sprains & injuries to tendons



VELSPAN™ (STRONG GRADUATED LONG STRETCH BANDAGE)

SIZE	CODE	QTY/BOX	BOXES / CASE
12cm x 4.5 m	770635	1	90

