



Datt Mediproducts Private Limited



# SOFTPAD™ (CDP)

## Sterile Combine Dressing Pad



### PROPERTIES

- Combine dressing of cotton pad with absorbing overwrap of non-woven material
- Individually packed in temper proof wraps
- Sterilized by Gamma irradiation/ ETO
- Non-adherent to the wound
- Highly absorbent cotton pad ensures slow strike-through even in case of heavy exudates

### COLOR

- Pure White

### MATERIAL

- Non Woven, Polyester



### INDICATION

- Wound dressing used to manage heavy exudates



### SIZES

7.5 cm x 7.5 cm

10 cm x 10 cm

10 cm x 15 cm

10 cm x 20 cm

