



Datt Mediproducts Private Limited



SOFTPAD™ (EYEPAD)

Eye Pad with Adhesive Borders



PROPERTIES

- Non-woven, Latex Free adhesive eye pad with adhesive borders
- Acrylic based Hypoallergenic adhesive
- Unique non-adherent absorbent pad coated with polyethylene mesh won't stick to the wound
- Gamma sterilization

COLOR

- White

MATERIAL

- Non Woven, Acrylic adhesive



INDICATION

- Eye injury care



SIZES

10 cm x 7.2 cm

