



Datt Mediproducts Private Limited



# VELDRESS<sup>®</sup>

## Non-Woven dressing with non-adherent pad



### PROPERTIES

- Convenient all in one dressing, reducing application time
- Non-Woven backing surface allows the wound to breath and water vapour to escape protecting skin from maceration
- Latex free hypoallergic adhesive
- Polyester wound contact layer prevents dressing from sticking to the wound & minimize pain during removal
- Absorbent non- adherent pad absorbs exudate & minimize fluid strike
- Leaves no traces of adhesive upon removal

### COLOR

- White

### MATERIAL

- Non- Woven, Acrylic Adhesive



### INDICATION

- Post-operative wounds
- Acute wounds (Abrasions & Laceration)
- Small or medium exudates level chronic wounds



### SIZES

5cm x 7.5cm  
10cm x 30cm

8cm x 15cm

10cm x 20cm

10cm x 25cm

