



Datt Mediproducts Private Limited

VELFIX[®] - EASY

Non-woven sterile IV dressing



PROPERTIES

- Sterile, Non-woven, IV dressing comes with two securing tapes
- Latex Free hypoallergenic adhesive used
- Highly porous & air permeable
- Protects injection site against dust, germs and other outside contaminants

COLOR

- White

MATERIAL

- Non-woven, acrylic adhesive



INDICATION

- I.V cannula fixation



SIZES

6 cm x 7 cm (S)

7 cm x 9 cm (M)

8.5 cm x 10 cm (L)

