



Datt Mediproducts Private Limited



VELFIX® - EASY+PAD

Non-woven dressing with Non-adherent pad



PROPERTIES

- Unique non-adherent absorbent pad coated with polyethylene mesh won't stick to the wound and absorbs small amounts of exudate, allowing normal wound healing process to occur with less pain and trauma
- Micropores in the non-woven material structure provides an adequate rate of moisture vapour permeability
- Covered with hypoallergenic acrylic adhesive
- ETO Sterile

COLOR

- White

MATERIAL

- Non-woven, acrylic adhesive



INDICATION

- Post operative wounds
- Acute wounds (Abrasions & Laceration)
- Small or medium exudates level chronic wounds



SIZES

5 cm x 7.5 cm
8 cm x 15 cm
10 cm x 35 cm

6 cm x 10 cm
10 cm x 20 cm

10 cm x 10 cm
10 cm x 25 cm

9 cm x 12.5 cm
10 cm x 30 cm

