



Datt Mediproducts Private Limited



VELFIX[®] - EDGE - S

Transparent I.V Dressings with border with silicone adhesive



PROPERTIES

- Soft silicone adhesive offers clean removal without causing pain
- Transparent Film to allow continuous observation of IV site
- Offers waterproof sterile barrier to external contaminants
- Notched design to discourage edge lift
- Comes with 2 securing tapes & a pre-printed documentation label
- Easy to handle with gloves allowing single handed application
- Specially designed for fragile & sensitive skin of neonates & elderly patients
- Ensures gentle yet effective securement
- Minimal stripping of skin cells & hair follicles upon removal
- Delivers constant adhesive strength
- Softer & more comfortable as compared to traditional adhesives

COLOR

- White with transparent PU film

MATERIAL

- Non-Woven Border
- Poly-urethane Transparent Film



INDICATION

- As Sterile I.V Dressing for cannula fixation



SIZES

5 cm x 6cm

6 cm x 7 cm

7cm x 9cm

