



Datt Mediproducts Private Limited



VELFIX[®] - EDGE

Transparent I.V Dressing with Borders



PROPERTIES

- Window frame design to allow continuous observation of IV site
- Offers waterproof sterile barrier to external contaminants
- Stabilization non-woven borders to maximize securement, breathability & wear time
- Notched design to discourage edge lift
- Comes with 2 securing tapes & a pre-printed documentation label
- Easy to handle with gloves allowing single handed application
- Unique coat technology with our Premium PU film offers provides MVTR of 2685 GSM/24 hrs
- EO Sterile

COLOR

- White with Transparent Window

MATERIAL

- Non-Woven Border
- Poly-urethane Transparent Film

INDICATION

- As Sterile I.V Dressing for cannula fixation



SIZES

5 cm x 6 cm

6 cm x 7 cm

7 cm x 9 cm

8.5 cm x 10 cm

10 cm x 12 cm

