



Datt Mediproducts Private Limited



VELFIX[®] - IV KIT

IV Cannulization Dressing Kit



PROPERTIES

- High utility dressing kit for peripheral lines designed to take care of insertion & dressing change of peripheral devices

CONTENT

- **Latex-free non-reusable tourniquet:** For finding veins singular usage, reducing cross- contamination due to carrier property of Velcro tourniquets and prevents skin allergy.
- **70% IPA Swab/CHG swab:** Is considered ideal for prepping the insertion site
- **Velfix Edge IV Dressing (7cm x 9cm):** Transparent IV dressing with Window Frame delivery design to allow continuous observation of the IV site while providing better seal around the catheter with stabilization non-woven boards
 - Notched design to effectively stabilize the catheter
 - Transparent film serves as a waterproof sterile barrier to external contaminants
 - Breathable transparent film provides high MVTR with an effective wear time of upto 7 days
 - One-Handed application & gentle removal
 - Made with a hypoallergenic, latex-free adhesive that is gentle to the skin yet holds catheter in place
- **Sterile Gauze Swab:** To control backflash of blood



INDICATION

- Sterile IV Dressing Kit for cannula fixation

