



Datt Mediproducts Private Limited



# VELFIX<sup>®</sup> - T+PAD

## Transparent film dressing with Non-adherent pad



### PROPERTIES

- Offers a transparent, waterproof film that serves as a bacterial barrier yet is breathable
- 30 µ thick PU film ensures increased breathability & softness
- Unique coat technology with patterned adhesive offers high MVTR of 2685 GSM/24 Hrs
- Latex free, non-adherent pad does not adhere to the wound bed, reducing pain upon removal
- Easy to mould around corners and is patient & clinician friendly
- Gamma Sterile

### COLOR

- Transparent

### MATERIAL

- Polyurethane film, Absorbent Pad, Acrylic adhesive



### INDICATION

- Post operative wounds
- Acute wounds (Abrasions & Laceration)
- Small or medium exudate level chronic wounds



### SIZES

|               |               |                |               |
|---------------|---------------|----------------|---------------|
| 5 cm x 7.5 cm | 6 cm x 10 cm  | 9 cm x 12.5 cm | 8 cm x 15 cm  |
| 10 cm x 10 cm | 10 cm x 20 cm | 10 cm x 25 cm  | 10 cm x 30 cm |
| 10 cm x 35 cm |               |                |               |

