



Datt Mediproducts Private Limited



SOFTLAPS[®]-Pw

Laprotomy Sponges



PROPERTIES

- Made out of Pre-washed 100% cotton gauze with no fluorescence
- Sponges are soft & have high absorbency
- Lint Free
- All edges are folded in and stitched with cotton threads, preventing the possibility of loose threads
- Sponges have dual- coloured radiopaque handles that remain visible when saturated with blood.
- Banded in steam compatible pack of 5

COLOR

- White with blue strips

MATERIAL

- 100% Cotton



INDICATION

- In OT to absorb blood



SIZES

25 x 25 cm x 8/12 ply

30 x 30 cm x 8/12 ply

40 x 25 cm x 8/12 ply

