



Datt Mediproducts Private Limited



SOFTLAPS®

Laparotomy Sponge



PROPERTIES

- Made of 100% Cotton Gauze Fabric
- Bleaching done by using Hydrogen Peroxide & de-mineralized water
- All four edges are sealed to prevent lint residue
- Available with or without X-ray detectable thread
- Tape stitched at the corner
- Available in sterile and non-sterile packaging
- Tied in bundles of 5's with red thread for easy counting

COLOR

- White

MATERIAL

- 100% Cotton



INDICATION

- In OT to absorb Blood



SIZES

20 cm x 20 cm x 4/6/8/12/16 ply
25 cm x 40 cm x 4/6/8/12/16 ply

25 cm x 25 cm x 4/6/8/12/16 ply
40 cm x 40 cm x 4/6/8/12/16 ply

30 cm x 30 cm x 4/6/8/12/16 ply
45 cm x 45 cm x 4/6/8/12/16 ply

