



Datt Mediproducts Private Limited



SOFTSWAB®

Sterile Gauze Swabs E.P



PROPERTIES

- Made of 100% Cotton Gauze Fabric
- Conform to EP- 17 Standards
- All swabs are with folded edges (American Fold), therefore no loose threads
- Available with or without X-ray detectable thread
- Bleaching done by using Hydrogen Peroxide and de-mineralized water
- ETO/ Gamma sterilization

COLOR

- White

MATERIAL

- 100% Cotton



INDICATION

- In OT to absorb blood
- Wound dressing



SIZES

5 cm x 5 cm x 8/12/16 ply 7.5 cm x 7.5 cm x 8/12/16 ply 10 cm x 10 cm x 8/12/16 ply

