



Datt Mediproducts Private Limited



VELSWAB®

Unsterile Gauze Swabs E.P



PROPERTIES

- Made of 100% Cotton Gauze Fabric
- Conforms to E.P-17 Standards
- All Swabs are with folded edges (American Fold), therefore no loose threads
- Available with or without X-ray detectable thread
- Bleaching done by using Hydrogen peroxide and de-mineralized water

COLOR

- White
- Green

MATERIAL

- 100% Cotton



INDICATION

- In OT to absorb blood
- Wound dressing



SIZES

5 cm x 5 cm x 8/12/16ply 7.5 cm x 7.5 cm x 8/12/16ply 10 cm x 10 cm x 8/12/16ply

Pack of 100's or Steam compatible pack of 2's/4's/5's with or without sterilization indicator or as desired by the Customer

