



Datt Mediproducts Private Limited



# Osteopop™

## Plaster of Paris Bandage

### PROPERTIES

- Premium quality of plaster of Paris with extra white finish
- Fast Drying
- Rich creamy Finish
- Early, easy & excellent mouldability and higher stability
- Interlocked weave of the cloth helps in plaster retention & provides high degree of accuracy
- Setting time is just 1.5-2 minutes
- Made from high purity gypsum, which gives rapid setting & durability
- Provide excellent finish & easy handling

### COLOR

- White


### MATERIAL

- Woven Cloth, Plaster of Paris



### INDICATION

- For immobilization to protect limbs with broken bones



### SIZES

10 cm x 2.7 m

15 cm x 2.7 m

