



Datt Mediproducts Private Limited



VELCAST™

Orthopaedic Fiberglass Casting Tape



PROPERTIES

- Excellent strength, rigidity & durable
- Good setting time & smooth moulding
- Excellent X-ray radiolucency
- Easy to apply and remove
- Lighter and stronger than the traditional POP bandage
- Water resistant
- Conforms to the area of body smoothly without wrinkles
- High quality knitted substrate and special resin protects against unexpected sudden impact
- Functional strength achieved in 20 minutes. Cost effective

COLOR

- White, Pink, Blue & Green (5 cm cast is available in polyester material with graphics)

MATERIAL

- Fiberglass, Polyester



INDICATION

- For Immobilization



SIZES

5 cm x 3.6 m

7.5 cm x 3.6 m

10 cm x 3.6 m

12.5 cm x 3.6 m

15 cm x 3.6 m

