



Datt Mediproducts Private Limited



Velnet®

100% cotton stockinette



PROPERTIES

- Protects the skin surface under compression
- 100% cotton for non-irritative comfort
- Easy application
- Lightweight and strong protection
- Conforms to joints, will not wrinkle or slip
- Durable and economical, may be washed for reuse
- Extra protection for frequent bandage changes
- Highly absorbent and well tolerated
- Steam autoclavable for sterile applications

COLOR

- White


MATERIAL

- 100% Cotton



INDICATION

- Protect skin from friction under plaster/ fibreglass bandage



SIZES

5 cm x 1 m
10 cm x 10 m

10 cm x 1 m
15 cm x 10 m

15 cm x 1 m

5 cm x 10 m

