



Datt Mediproducts Private Limited



SILVEL™

Nanocrystalline Silver coated Antimicrobial Barrier Dressing



PROPERTIES

- Three layered dressing. Wound contact layer is made of 100% Cotton, middle layer is made of Viscose and outer layer is made of polurethane net
- Sustained release of Ag⁰ is 294 ppm/ 24 hours
- Ag⁰ content in Silvel is 438 mg/100 cm²
- Provides long lasting anti-microbial barrier (approx. 5 days)
- Effective against resistance bacteria like MRSA/VRE/Pseudomonas
- Maintains moist wound environment optimal for fast healing
- Silvel needs to be moistened only once prior to use, repeated moistening is not required
- Non adherent in nature
- Gamma Sterile

COLOR

- Amber

MATERIAL

- Nanocrystalline Silver coated non adherent cotton fabric, Absorbent Pad



INDICATION

- Partial thickness burns
- Pressure ulcer
- Diabetic foot ulcer
- Infected wounds
- Surgical wounds
- Venous leg ulcer



SIZES

5 cm x 5 cm	10 cm x 10 cm	10 cm x 20 cm	15 cm x 15 cm
20 cm x 20 cm	20 cm x 40 cm	40 cm x 40 cm	10 cm x 120 cm

