



Datt Mediproducts Private Limited



# MINIPORE™

## Hypoallergenic & Latex free Surgical Paper Tape



### PROPERTIES

- High quality paper with gentle adhesive
- Hypoallergenic and latex free
- White, lightweight, non-woven paper tape
- Excellent grip
- Porous structure allows skin to breathe
- Conforms well on skin
- Easily hand tearable
- No traces of adhesive remains on the skin after removal

### COLOR

- White

### MATERIAL

- Non Woven, Acrylic Adhesive



### INDICATION

Perfect for securing of:

- Dressing
- Tubes
- Catheters, etc.



### SIZES

1/2" x 9.1/5 m

1" x 9.1/5 m

2" x 9.1/5 m

3" x 9.1/5 m

