



Datt Mediproducts Private Limited



VELPORE™ -S

Hypoallergenic tape with soft Silicone Adhesive



PROPERTIES

- Offers extra securement of dressings & tubes.
- Constant adhesive level throughout the wear time.
- Gentle to skin, Hypoallergenic
- Breathable, Non-woven paper
- Water-resistant
- Can be Repositioned
- Latex free, Light weight
- Easily tearable by hand
- Clean removal without causing pain

MATERIAL

- Non -Woven, Soft Silicone Adhesive

BENEFITS OF SILICONE ADHESIVE

- Softer and more comfortable
- Lower surface tension
- Works differently with skin and hair
- Naturally conforms to skin profile
- Doesn't stick to gloves

APPLICATION

- I.V. Tubing
- Probes
- Cannula
- Ideal for sensitive skin, Neonates, elderly patient

INDICATION

- Paediatrics & Neonatologist, Dialysis centre, ICU, Medical Surgery, Emergency

SIZES

Product Name	Code	Size	Qty/Box	Boxes/Case
Velpore-S (Soft Silicone Adhesive Tape)	770863	1 inch x 5m	12 Rolls	6 Boxes
Velpore-S (Soft Silicone Adhesive Tape)	770864	2 inch x 5m	6 Rolls	6 Boxes
Velpore-S (Soft Silicone Adhesive Tape in individual box pack)	770865	1 inch x 1.4m	50 Rolls	4 Boxes
Velpore-S (Soft Silicone Adhesive Tape in individual box pack)	770866	1 inch x 1.4m	25 Rolls	4 Boxes

