



Datt Mediproducts Private Limited



VELSILK™

Hypoallergenic & Latex free Silk-like Tape



PROPERTIES

- Latex-free, hypoallergenic silk-like tape
- Adheres well with good initial stick, High strength
- Water resistant

COLOR

- White

MATERIAL

- Silk-like cloth, Acrylic adhesive



INDICATION

- Fixation of dressings on sensitive skin
- Securing of :
- I.V Tubing
 - Catheters
 - Probes
 - Cannulae
 - Stabilizing finger splints



SIZES

1/2" x 9.14 m

1" x 9.14 m

2" x 9.14 m

3" x 9.14 m

SINGLE USE

1" x 1 m

2" x 1 m

