




Datt Mediproducts Private Limited

Velsoft® FR-T

Transparent Fixing roll for dressing fixation & skin protection



PROPERTIES

- Light Weight, Fully transparent, skin friendly
- Latex free Hypoallergenic Adhesive
- Easy to Remove
- Conforming & Durable
- Water resistant property, allows patients to shower
- Barrier to external contaminants
- Allows skin to breath thus reducing risk of wound infection.
- Leaves no traces of adhesive on the skin upon removal

COLOR

- Transparent

MATERIAL

- PU film , Acrylic adhesive



- Provide waterproof wide-area fixation of wound dressings, also on joints and extremities



SIZES

5 cm x 2 m

5 cm x 10 m

10 cm x 2 m

10 cm x 10 m

15cm x 2 m

15cm x 10 m

