

EZ-CLEAN™ RECOVERY HEART PILLOW

A WATERPROOF THERAPEUTIC TOOL & PERSONALIZED GET-WELL GIFT

FREE TRIAL
30 DAYS



A GERM-RESISTANT BRAND HOSPITALS TRUST

Post-Op Support
For Patients after Surgery
on Heart, Lung, Breast,
Stomach, Etc.

Vomit Support
For Expectant Mom and Ones
with Stomach Flu, Acid Reflux
and Heartburn Concerns

Spill Support
For Families with
Children that Have
Spills and Messes





EZ-CLEAN RECOVERY HEART PILLOWS CREATED TO REDUCE INFECTION RISK



EZ-Clean Recovery Heart Pillows is a trending brand developed in 2016 after Baystate Heart Hospital brought us a common problem in desperate need of a solution. The recovering heart surgery patients are exposed to infection, even the post-operative complication when they hug a cloth-surfaced heart pillow to protect the incision site as a necessary sternal precaution or a coughing exercise.



Cloth-Surfaced Pillows are Germ's Residence

Messes can easily stick to the rough surface of a cloth pillow. Any liquid, including vomit, blood, drinks, or water can seep between the threads, leak inside, and soak the fiberfill. Germs on the uneven threads and the inner fiberfill are not removable, but can easily spread to the patient users. Shown to the left is a recovering cardiac survivor who is embracing a porous cloth-surfaced heart pillow against his chest, just an inch from the incision site.

Vinyl-Covered Pillows are Germ-Resistant

Alongside the hospital, we innovated the EZ-Clean Pillow with a removable waterproof cover that wraps the inner pillow. Equipped with the EZ-Clean Sanitary Solution as illustrated below, our pillow is much safer. It resists germs, stops the spreads of germs, and ensures patient safety effectively.



INFECTION CONTROL MADE EASY AND EFFECTIVE WITH OUR EZ-CLEAN SANITARY SOLUTION



OUTER COVER

- Made of specialty vinyl
- Smooth surface
- Waterproof, wipeable, and washable
- Removable
- Durable and ink friendly
- Firm to support, comfy to touch



INNER PILLOW

- Stays clean and dry
- Filled with micro-beads
- Beads move to contour to the user's body upon pressed offering custom support effect



SAFE OPENING

- Waterproof protective flap
- Controlled by a hidden zipper
- Cloth covered zipper ends
- Unintentional access to inside prevented



WATERPROOF

- No vomit, blood, mucus, other liquid, mess, or germs can break through the cover
- The waterproof cover keeps the inner pillow dry and clean making the entire pillow last longer.




WIPEABLE

- To clean the germ-resistant vinyl cover, wipe the smooth surface with tissue or cloth
- Continue using the pillow safely after spraying a disinfectant



WASHABLE

- To sanitize the pillow, wash and dry the removable cover
- Wash by hand or machine
- Hang dry or use a drier
- The cover can last at least 15 washing cycles



AN EZ-CLEAN THERAPEUTIC TOOL AND PERSONALIZED GET-WELL GIFT



EZ-Clean Heart Pillows are waterproof therapeutic tools (post-op support, vomit support and spill support), and personalized get-well gifts for recovering patients who had a surgery on heart, lungs, breasts, liver, stomach, intestine, kidneys, as well as supports designed for moms in and following labor. EZ-Clean Pillows are also ideal for those who are prone to vomiting such as pregnant women and anyone else with the stomach flu, acid reflux or heartburn, and those who prefer spill-proof, hygienic, hypo-allergenic comfortable cushion. Visit www.bestplush.com/users-guide for a detailed Guide.



Splint Incision/Sternal Precaution

After the surgery when the bone heals, patients must protect it from the excess force and limit movement of the chest and upper body. Hugging our pillow as guided can splint their external incision and conduct necessary movements safely.



Coughing and Deep-Breathing Exercise

Recovering from a cardiac, chest or esophageal surgery, patients often avoid coughing or deep breathing which causes pain. This can build up mucus and phlegm in the lungs, slowing down blood circulation and oxygen moving into the lungs, which can cause pneumonia and other lung ailments. Coughing and deep breathing exercises aid in prevention.

Patients shall place the Heart Pillow with the V-cut towards the chin, wrap their arms around the pillow, hold it tightly, cough and breathe deeply as guided. Proper pressure placed on the incision site will help control and reduce jolting and pain.



Maintain Positions

- Supine position
- Upright sitting position
- Fowler's position
- Healthy sitting position on bed and in chair



Protect from Vomit, Spill, and Liquid



Personalized Encouragement Gift That Lasts Long

Families, friends, doctors, and nurses can sign the smooth, germ-resistant cover. The bright, get-well messages provide strong emotional support and encouragement to the patients, particularly for those who are experiencing a prolonged health event. The personalized heart pillow helps them stay positive, aiding in a quick recovery and healthy life. After recovery, the therapeutic tool will become an unforgettable keepsake and easy-clean cushion for home use.



COMPARE BEFORE BUYING

FREE TRIAL
30 DAYS

LET'S COMPARE	OUR HEART PILLOWS	OTHER HEART PILLOWS
INCLUDES	2 pieces: outer cover and inner pillow	1 piece: pillow only
OUTER COVER	Vinyl outer cover, removable	No outer cover
SURFACE MATERIAL	<ol style="list-style-type: none"> 1. Outer cover made of vinyl 2. Waterproof 3. Smooth 4. Durable 5. Firm yet comfortable, more supportive 6. Removable 	<ol style="list-style-type: none"> 1. Surface made of cloth 2. Porous 3. Rough, uneven 4. Less durable 5. Soft, comfortable, less supportive 6. Not removable
FILLING MATERIAL	Microbeads	Fiberfill
WATERPROOF	Waterproof cover prevents leaks, keeping inner pillow dry and clean	Fluids and germs penetrate the fabric outer layer affecting the inner fill
CLEAN AND REUSE	<ol style="list-style-type: none"> 1. Remove mess from the smooth cover by simple wipes and spray with a disinfectant 2. It is safe to continue using it immediately 3. Backup pillow is not necessary 	<ol style="list-style-type: none"> 1. Rough absorbent fabric surface and inner fiberfill prevent cleaning and disinfection 2. It is not safe to continue using a messy pillow 3. Backup pillows are necessary
WASH AND DRY	<ol style="list-style-type: none"> 1. Wash/dry the removable cover only 2. Wash by machine or hand 3. Hang dry or use a drier 4. Stay in tact at least for 15 washes 	<ol style="list-style-type: none"> 1. Wash/dry the entire pillow 2. Machine wash only 3. Machine dry only 4. Fiberfill begins hardening by the third wash
SUPPORT	Micro-beads move to contour to the user's body when hugged so the user can adjust the hugging force to customize supportive effect	Fiberfill does not contour to the user's body
PERSONALIZE	<ol style="list-style-type: none"> 1. It is easy to print logo, sign names, and write messages on the vinyl cover 2. Ink dyes surface only and shows bright 3. Ink stays clear for years, even through cleaning and washing 	<ol style="list-style-type: none"> 1. It is difficult to print or write on the soft rough cloth surface 2. Ink sinks into the weaves, and shows fuzzy 3. Ink fades quickly, particularly in washing or cleaning

Our EZ-Clean Recovery Heart Pillows were developed alongside a heart hospital to replace the cloth-surfaced heart pillows with infection prevention as the top priority. Compared to the cloth pillows, our waterproof, wipeable and washable pillows are much easier to clean, safer to use, more supportive and durable. That is why our pillows have been trending and rapidly replacing the cloth-surfaced heart pillows in many major hospitals since Covid-19 occurred. After the pandemic, patients and non-patient consumers will surely further value germ-resistant products like our heart pillows that hospitals trust.

If you are looking for an easy-clean post-operative tool, comfortable therapeutic aid, and personalized get-well gift idea for your customers, contact us for a free sample or order an Introductory Kit for a 30-day free trial.



- Large: 17.5" x 16.5"
- Medium: 12.5" x 12.0"
- Price at \$11.99 and up
- Minimum order for 10 pillows
- Get-Well-Soon optional print
- Custom order welcomed
- Ready to Ship from NJ 08036

USES	USERS	SIZE
Post-Op Support	Surgery on heart, lungs, breasts, liver, stomach, intestines, kidneys and for moms in and following labor	Large & Medium
	Breast cancer surgery; mastectomy	Medium
Vomit Support	Pregnant women and those with stomach flu, acid reflux or heartburn	Large & Medium
Spill Support	Those who prefer spill-proof, hygienic, hypo-allergic comfortable cushion	Large & Medium